silasoft® Normal

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silasoft® Normal

Instructions for use

Precision impression material, C-silicone based, light bodied consistency

Fig. 1

1. Preparation of the impression

Depending on the impression technique select an appropriate tray. For optimal adhesion of the impression to the tray apply a thin layer of **sili Adhesive Spray** or **sili Adhesive liquid** onto the tray and let dry for 60 sec (Fig. 1). Take preliminary impression with **silaplast FUTUR**. To facilitate the flowing out of excessive light body material from the impression cutting of escape vents and removal of interdental septi in the impression with the use of the special instrument **Deta-Cut** is recommended.

Fig. 2

2. Mixing

Extrude **silasoft® Normal** and **cat® universal FUTUR** in equal lengths and thickness onto mixing pad (Fig. 2). Using catalyst liquid, the mixing ratio is 1 cm of **silasoft® N** to 1 drop of catalyst. Mix according to dosing ratio for approx. 30 sec. with spatula until homogeneous uniform colouring is achieved. Load preliminary impression. A normal syringe or disposal syringe can be used for bubble free application (Fig. 3 + 4).

Fig. 3

3. Disinfection

After removal rinse impression under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

4. Model casting

Casting between 30 minutes and 24 hours after removal of the impression. Recommended materials: Dental plasters of class III and IV - also commercially available model acrylics.

Fig. 4

5. Galvanisation

Impressions can be electroplated with commercially available copper and silver plating baths.

6. Cleaning of trays

Remove set material with a blunt instrument. **sili Spray** or **sili liquid** film can be easily removed with **sili Spray Solvent** (free of FCHC). Work in well ventilated rooms. Clean und disinfect as usual.

Important working hints


- The good flowability of **silasoft® N** guarantees an outstanding detail reproduction. Due to this very light-flowing consistency however, **silasoft® N** tends to a partial segregation in rare cases. Please homogenize the eventually segregated **silasoft® N** with a spatula on a mixing pad before mixing with the catalyst.
- Increased temperatures accelerate, decreased temperatures retard the setting time.
- Catalyst bottles or tubes should be closed tightly immediately after use.
- Don't leave any residual material in the patient's mouth.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Avoid contact of catalyst with eyes, may cause irritation. If necessary flush eyes with plenty of water and seek immediate medical care.
- Please follow the instructions of the safety data sheet!

Avoid contact with eyes and skin!
Wear suitable eye protection and gloves!

Other information:
Silicone impression materials are proven a million times over, unwanted effects cannot be expected in proper use. However, immune reactions, e.g. allergies or irritations cannot be excluded basically. In the case of doubt, we recommend to get tested the allergenic potential before the application of the material.

Only for dental use by qualified personnel.

silasoft® N / cat® universal FUTUR:

Indications for use:

- Corrective ("Wash") material for double impression
- Reline impression

Technical Data:

DIN EN ISO 4823 - Type 3

- **Mixing ratio:**
silasoft® N / cat® universal FUTUR
1 cm silasoft® N = 0.50 g / 1 cm cat® = 0.05 g
silasoft® N / silasoft® cat®
1 cm silasoft® N = 0.50 g / 1 drop cat® = 0.02 g
- **Colour code:**
silasoft® N: rose-opaque
cat® universal FUTUR: blue
silasoft® cat®: light blue
- **Mixing time:** approx. 30 sec.
- **Working time:** 1 min. 15 sec.*
- **Setting time:** approx. 3 min. 45 sec.*
- **Time in mouth:** approx. 2 min. 30 sec.
- **Strain in compression:** approx. 11.2 %
- **Recovery from deformation:** approx. 98.4 %
- **Linear dimensional change:** approx. 0.5 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:** 15 °C / 59 °F to 25 °C / 77 °F

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering information:

- silasoft® Normal Standard packing** 02004
160 ml tube
- Tube 4-pack** 02005
4 x 160 ml tubes
- Standard bottle-packing** 02320
160 ml bottle
- Bottle 4-pack** 02321
4 x 160 ml bottles
- cat® universal FUTUR** 02015
5 x 35 ml paste
- silasoft® cat®** 02017
5 x 10 ml bottles

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 07243/510-0 · Fax: 07243/510-100
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personell).

Made in Germany


11/2017

silasoft® Normal

Precision impression material, silicone based, condensation curing

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