onetime perfect wash

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onetime perfect wash

GB
Instructions for use

Precision impression silicone, addition curing, light bodied consistency- fast setting

Fig. 1

1. Tray preparation

Select a tray type suited for desired application (standard-size or custom, in metal or hard plastics). To avoid teeth contact with the tray, place stops from wax or light-curing resin into the tray at non-relevant areas of the impression. For optimal adhesion of impression material, we recommend applying a thin layer of **trayloc® A** to tray (Fig. 1). Let dry for approx. 3 minutes.


Fig. 2

2. Mixing and dispensing

Insert cartridge into Automix2-System mixing gun (Fig. 2). Remove and dispose of closure cap. For exact flow control extrude slightly material until uniform flow from both orifices is achieved (Fig. 3). Mount the mixing cannula by using the guides on the cannula and the cartridge and lock in place by rotating in the opposite direction (Fig. 4). Feed material applying even pressure. Prior to use, dispense a small amount of material and visually check that the base and catalyst paste are mixed omogeneously (Fig. 5). Only then dosing can be performed individually. After each use leave the used mixing cannula in place to seal the cartridge.


Fig. 3

3. Using the sandwich technique

First apply the mixed **onetime perfect putty** onto the impression tray, giving the teeth area a distinctly trough-like shape (Fig. 6). Then apply a layer of **onetime perfect wash** into the shaped depression (Fig. 7). For applications in the groove area use intra-oral-tip or disposable syringe (Fig. 8). Insert the filled tray into patient's mouth and position it within 3 – 4 seconds, under light pressure. Then seat the impression without pressure, until **onetime perfect** has set.


Fig. 4

4. Disinfecting

After removing tray from the mouth, rinse impression under running, lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde or other usual solutions for 15 minutes.


Fig. 5

5. Pouring the model

We recommend waiting at least 30 minutes before pouring the model. The impression will retain its dimensional stability for up to 14 days. All class III and IV dental plasters and standard modelling resins can be used.

Fig. 6

6. Cleaning the tray

After setting, material can be removed with a blunt instrument. Remaining film of **trayloc® A** is dissolved by immersion in acetone or isopropanol (be sure to use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.


Fig. 7

important working hints

- Substances such as retraction liquids can adversely affect the setting reaction and should therefore be thoroughly rinsed off before taking the impression.
- Do not combine with condensation curing silicones.
- Contact with latex examination gloves or other latex material may impair the setting reaction of **onetime perfect**.
- Be sure to remove all traces of impression material from patient's mouth.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- We recommend use of **onetime perfect putty** as tray material in the double-mix technique.
- Please follow the instructions on the safety data sheet!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, **undesired effects** are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Only for dental use by qualified personnel.



Fig. 8

Indications for use:

- Impressions with double-mix and sandwich technique
- Crown- / bridge work
- Inlays, onlays, partial crowns
- Dowel crowns, veneers

Technical data:

- DIN EN ISO 4823 - Type 3
- **Mixed volume:** 50 ml cartridge
 - **Mixing ratio:** 1:1
 - **Product colour:** base: violet catalyst: purple
 - **Mixing time:** omitted (Automix2-System)
 - **Working time:** approx. 1 min. 30 sec.*
 - **Setting time:** approx. 3 min. 30 sec.*
 - **Time in the mouth:** approx. 2 min.
 - **Strain in compression:** approx. 6,5 %
 - **Recovery from deformation:** approx. 99,8 %
 - **Linear dimensional change:** approx. 0,2 %
 - **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity.
 - **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

Ordering information:

- onetime perfect Introkit 02619**
2 jars of putty, ea. 250 ml
2 cartridges of wash, ea. 50 ml
2 measuring scoops
12 mixing cannulas, pink
12 intra-oral-tips, white
10 disposable impression syringes, 1 flowchart

- onetime perfect wash Standard-Pack 03289**
2 cartridges, ea. 50 ml
12 mixing cannulas, pink
12 intra-oral-tips, white

- Eco-Pack 02152**
8 cartridges, ea. 50 ml
48 mixing cannulas, pink
48 intra-oral-tips, white

- onetime perfect putty Standard-Pack 03288**
2 jars, ea. 250 ml
2 measuring scoops

- Eco-Pack 02155**
8 jars, ea. 250 ml
4 measuring scoops

- trayloc® A 03098**
17 ml bottle with brush
15 ml spray bottle **03392**

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