# QUBE® /Plus/Premium

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D

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Original  
Instructions

E

We are pleased that you decided to buy a highly developed piece of equipment from SCHICK and would like to wish you every success when working with your new **QUBE** control unit.

We wrote these instructions to enable you to get accustomed to your new piece of equipment and to provide you with the correct operating and maintenance instructions.

Index:	page
1. Safety information.....	2
Caution!	
2. Range of applications.....	2
3. Summary - List of contents.....	3
4. Accessories/options.....	4
5. Initial start-up and installation switching on and off.....	5-6
6. Control satellite.....	7
7. Tool change on motorhandpiece.....	8
8. Operation of touch-sensitive surfaces..	9
9. Basic and additional functions.....	9-11
10. Menu structure and basic settings.....	11-14
11. Speed control functions.....	15
12. Connecting and controlling a second.. motor handpiece	15
13. Cleaning and Maintenance/..... removing the chuck	16-17
14. Possible error messages.....	18
Reset after error message	18
15. Suggestions for wax knife.....	18
16. Technical Data.....	19
17. Declaration of conformity.....	20

## 1. Safety information

- 1.1 Ascertain that your mains supply coincides with the data on the rating plate.
- 1.2 **QUBE**- units are not suitable for the following applications:
  - in areas where there is a risk of explosion
  - for medical applications
  - for working on moist materials
- 1.3 Ensure that all regulatory requirements are observed during use
  - always wear protective glasses
  - provide enough light at the working place
  - use dust suction
- 1.4 Under no circumstances should the motor handpiece be cleaned with compressed air.
- 1.5 Before putting the handpiece down, always insert a rotary instrument or the pin, supplied with the unit, into chuck.

### Caution!

- When using rotary instruments, do not exceed the maximum speeds laid down by their manufacturer.
- When operation with left-hand rotation, the collet may loosen when there is a high load.
- Repairs and other technical procedures must only be carried out by suitable qualified personnel, authorized by SCHICK.
- SCHICK do not guarantee the **QUBE** unit should it not have been used in accordance with the instructions.

## 2. Range of applications

The **QUBE** unit is designed for universal use in dental laboratories when trimming crowns and bridges, respectively acrylic and light chrome cobalt dentures. With its speed range from 200 - 60.000 rpm the SCHICK **QUBE** - unit allows to work all dental materials.

The motorhandpiece is ergonomical, small shaped and runs smoothly, thus eliminating user-fatigue and making it economical to use over long periods. The progressive processor controller constantly monitors the unit for overloading. A second thermal fuse, which is independent of the computer, provides further protection. Thus the **QUBE** - unit offers a maximum of safety.

### Conditions of environment:

- interior 5° - 40° C.
- up to 2,000 meter over sea level


Categorie of overvoltage: II

Grade of pollution: 2

WEEE-Reg.-Nr. DE 78620387

### 3. Summary **QUBE** - List of contents

knee control:

Art.no.

<b>QUBE</b> knee controller	(1a)	9016
<b>QUBE</b> Plus		9036
<b>QUBE</b> Premium		9056
motorhandpiece with cable	(2a)	9001
motorhandpiece LV w. cable	(2b)	9002
control satellite	(3)	9020
handpiece rack with two keys for changing chuck	(4)	9127
power supply	(5a)	9102
velcro	(5b)	9130
mains cable	(6)	2160
suspension strip	(7)	9103
screws - 2 pieces-	(8)	3170
cable satellite 2m	(9)	9124
signal transmitter for suction	(10)	9060

foot control:

Art.no.

<b>QUBE</b> foot controller	(1b)	9017
<b>QUBE</b> Plus		9037
<b>QUBE</b> Premium		9057
motorhandpiece with cable	(2a)	9001
motorhandpiece LV w. cable	(2b)	9002
control satellite	(3)	9020
handpiece rack with two keys for changing chuck	(4)	9127
power supply	(5a)	9102
velcro	(5b)	9130
mains cable	(6)	2160
cable satellite 2m	(9)	9124
signal transmitter for suction	(10)	9060

table top control:

<b>QUBE</b> control unit table top model	(1c)	9018
<b>QUBE</b>		9038
<b>QUBE</b>		9058
motorhandpiece with cable	(2a)	9001
motorhandpiece LV w. cable	(2b)	9002
control satellite	(3)	9020
handpiece rack with two keys for changing chuck	(4)	9127
power supply	(5a)	9102
velcro	(5b)	9130
mains cable	(6)	2160
suspension strip	(7)	9103
screws - 2 pieces -	(8)	3170
cable satellite 2m	(9)	9124
signal transmitter for suction	(10)	9060

## 4. Accessories/options

chuck Ø 2,35mm		
foot switch (for table top model only)	(11)	9142 / 6642 (LV)
dynamic foot-switch (for table top model only)	(12)	6370/2 9440
chuck key		
engineers wrench		4115
motor cable complete		9188 / 6223 (LV)
cable satellite 0,5m		9190 / 9225 (LV)
cable satellite 1,5m		9124/1 9124/2

riveting hammer(1)  
Art.no. 1850/2


Q Motor with ISO coupling and cable type 4  
art.-no. 9002/04


## 5. Initial start-up and installation

### 5.1 Installation of the suspension strip for knee control unit or work bench control unit

1.

2.


- 1 The suspension strip (7) is supplied with a template which acts as an aid for correct positioning of the **QUBE** knee control unit. To do this, the mounting bracket is positioned on the work bench as shown in figure 1, and secured using the screws supplied (8). After installation, the template can simply be detached from the bracket.
2. For installation of the knee control unit, the cavity at the rear of the unit is placed on the bracket and pushed back until it clicks into place (Fig.2).

### 5.2 Fixing the power supply with the velcro

The delivered velcro can be used to attach the power supply under the workbench (Fig.2).

1.

2.


1. Remove one of the protective foils of the velcro and attach it to the power supply.
2. Remove second protective foil and stick the power supply to a clean and even surface. Please avoid permanent load to the power supply by the cables.

### 5.3 Initial start-up

Connect the network cable (6) to the power supply (5a). Connect the motor hand piece to the 'motor hand piece 1' socket. Connect the controls satellite (3) to the control unit using the cable (9). Plug the network adaptor into the power supply socket on the control unit. To switch on the device, activate the operating component (knee pad, foot lever, etc.) once.

If the worktop version is being used, the device is switched on by pushing the MENÜ touchfield.

### 5.4 Switching on and off


There is no main switch for the QUBE control unit. If the hand piece is not used for more than 60 minutes, it switches off automatically. Alternatively, the device can be switched off manually by simultaneously depressing the touch-field ESC + for longer than 2 seconds. By simply touching the control device (knee pad /foot pedal, MENÜ button on the worktop version), the device is switched on again.

 In order to prevent unnecessary waiting times, when the wax knife is switched on, the control mechanism is not automatically switched off.

 Basically, firstly attach the handpiece, the controls satellite and the electric wax knife (if you have one) before connecting the mains cable and switching the control mechanism on.

### 5.5 Connection of the signal transmitter for automatic switching on of a workplace suction device

In order to operate the workplace system together with a workplace suction system with automatic function, the signal transmitter for the suction signal must be connected. This is, as shown in the illustration, inserted in the suction signal jack on the control unit. At the opposite end, in the appropriate connector socket on the workplace suction device.


When operation QUBE together with a KaVo SMARTair suction or a Zubler suction AP 501, the optional cable Art.No.:9229 can be used as signal connection for controlling the suction.

## 6. Control satellite

### 6.1 The touch pad

On the top of the control satellite there are five touch-sensitive keys (hereafter referred to as 'touch-field') and a touch sensitive control dial (hereafter referred to as 'wheel').

### 6.2. Operation of touch-sensitive surfaces (touch-fields)

Operation of the controls satellite is carried out by a simple touch of the touch fields. Within this, there is a distinction between a single touch, holding for 2 seconds and a rotational motion of the wheel. Each action can be found in the appropriate section of this operating manual.

basic operation:


a circular motion on the wheel gives maximum torque, or the speed setting can be adjusted via the menu. Activating the touch-field **Menu** for two seconds opens the settings menu. One touch enables you to skip to the next menu point. Simply touching the four outer soft keys enables operation of the additional functions, activation of the cruise control or selection of the four program Memory within the basic display.


In the **QUBE** table top version, the hand piece is started by a single touch of the touch-field **Menu**

### 6.3 LCD- Information display

The symbols

### 7. Tool change on motorhandpiece


By turning the handle on the hand piece, the chuck can be opened or closed

 Only carry out a tool change when the motor is switched off! To optimise the precision and durability of the chuck, the tool or the factory-supplied protective pin must always be in a fixed position, even when not in use.

 Always insert tools as far as possible into the chuck in order to achieve maximum holding power!

## 8. Operation of touch-sensitive surfaces


There is continuously variable speed control of the handpiece up to the maximum preset motor speed using the respective control unit (see diagram below). Exception: using the footswitch that is attachable to the bench control unit, the handpiece is regulated directly to the preset speed when activated. In relation to this, please see the instructions given in point 9.

## 9. Basic and additional functions, time and thermometer

### 9.1 Pre-selection of maximum speed

1.

1. Adjust the maximum speed by making a rotational motion on the wheel.

For using the speedrange between 200 and 1.000 rpm, the maximum speed has to be limited to 1.000 rpm at the control satellite

Please always observe the maximum permitted speed for your tools!

Flashing of the rpm symbol illustrates the pre-selected rpm when the motor is not running. If the handpiece is started using the knee pad or foot lever, the display changes to the current rpm and the rpm symbol is displayed permanently.

When operating the motor hand piece without the controls satellite, the maximum speed is restricted to 30,000 rpm.


## 9.2 Additional functions, time and thermometer

Using the controls satellite, the following additional functions are available; time, date, stopwatch, countdown, room temperature.

After initial set-up of the unit, the time appears in the LCD display in the bottom left-hand corner, and parallel to it, the time symbol. These additional functions can all be called up one after the other using the 'additional function' touch field.

### 9.2 Time


1.

1. Basic time display

### 9.3 Date display


1.


1.  Press to call up date

### 9.4 Stopwatch


1.

2.


3.

1.  Press to call up stopwatch
2.  Press to start or stop stopwatch
3. Reset to 0, press  for 2 seconds


## 9.5 Countdown

1. Press to call up countdown
2. Enter the required time using the wheel
3. Press to start or stop countdown

Following expiry of the time entered, there is an acoustic signal

## 9.6 Room temperature display

1. Press to call up the current room temperature, as measured on the control satellite



One more touch of the time signal activates a return to the time display.

## 10. Menu structure/basic settings, time and date settings together with user-specific programming

To activate/deactivate the wax knife, adjust the date and time, adjust the direction rotation, and to save individual performance and speed programmes, you need to access the programming level. From each menu you can go to the basic display by using the touch-field . The settings are stored.


### 10.1 Switching on the wax knife

(in this respect, please note the recommendations on page 18 and the separate instruction sheet for the wax knife)

1. **Menü** Press menu for 2 seconds to switch to programming level
2. To activate the wax knife, set a value of 0 by turning the wheel  
The heating can be adjusted within the range of 1 to 20

### 10.2 Switching from right or left rotation



1. **Menü** Press menu to adjust from right or left rotation
2. Turn the wheel to select the desired direction of rotation of the motor handpiece


The chuck can become detached when rotating to the left under high strain

### 10.3 Torque limitation


1. **Menü** Press menu to access the torque limitation function
2. Adjust the maximum torque by making a rotational motion on the wheel

### 10.4 Selection of Motor 1 and Motor 2 (Only with **QUBE** Premium. For further information see point 12, page 15)


1. **Menü** Press menu to arrive at the motor selection function
2. Select motor 1 or Motor 2 by rotating the wheel

## 10.5 User-specific programmes

1. **Menu** Press menu to switch to programme selection
2. Set the programme to be changed by making a rotational motion on the wheel
3. **Menu** Press menu and adjust maximum speed using the wheel
4. **Menu** Press menu and adjust maximum torque using the wheel
5. **Menu** Press menu and select Motor 1 or Motor 2 with the wheel (only with Premium)
6. By pressing **Menu** again, return to menu programme selection. Here, you can now configure a further programme by using the wheel, or by pressing **ESC**, leave the programming level

## 10.6 Setting the date


1. **Menu** Press menu to switch to date settings
2. Adjust the date format by making a rotational motion on the wheel
3. **Menu** Press menu to adjust year and change with the wheel

4. **Menu** Press menu to adjust month and change with the wheel
5. **Menu** Press menu to adjust day and change with the wheel and then confirm with menu

### 10.7 Time setting

1. **Menu** Press menu to switch to time setting
2. Adjust the time format (12h/24h) by making a rotational motion on the wheel
3. **Menu** Press menu to adjust seconds and change with the wheel
4. **Menu** Press menu to adjust minutes and change with the wheel
5. **Menu** Press menu to adjust hours and change with the wheel and then confirm with menu

### 10.8 Operating hours counter

1. **Menu** Press menu to display the operating hours counter

After pressing **Menu** once more, the display returns to the basic display

Due to the break in cycle at Schick, up to 24 hours can be displayed.

It is possible to leave the programming level at any time by touching the touch field **ESC**. Changes already made are not lost as a result.

## 11. Speed control functions

The motor handpiece can be operated using the speed control function without having to continuously press the control unit.

Two different speed controls can be used via the controls satellites.

1. Simple speed control: speed is maintained without pressing the control lever further as soon as the speed has been held constant for longer than 2 seconds
2. Speed control with touch function: The motor handpiece can be run at the preset speed by a single touch of the operating lever.



Please always observe the maximum permitted speed for your tools!

To activate the speed control function, proceed as follows:

### 11.2 Simple speed control

1.



1.  Press to activate the simple speed control function.

The speed control is now active. If a selected speed is held for more than 2 seconds, this speed is saved to the memory and the operating unit (knee pad etc) is no longer needed. To stop, touch the operating unit once

### 11.3 Speed control with touch function

1.

1.  Press to activate speed control with touch function.

Speed control with touch function is now active. Once the control unit is pressed, the handpiece starts at the preset speed on the control satellite.

## 12. Connecting and controlling a second motor handpiece (only with **QUBE** Premium)

The **QUBE** Premium control unit enables a second motor handpiece to be connected. For this purpose, a second **QUBE** hand piece or a hand piece from the Profi, Junior or C Master ranges that can be connected using an adaptor cable (obtainable from Schick), would be suitable.


The handpieces can only be operated individually and not in parallel!

## 13. Cleaning and Maintenance / Removing the chuck

The **QUBE** motor handpiece is designed for maximum durability and therefore, the chuck and the tip should be removed and cleaned from time to time.

### 13.1 Removal of chuck

1. Remove tool, unscrew tip of handpiece and remove dustcap
2. Place the chuck key (Art. no. 4115) on the open chuck, and with an open spanner, (Art.no. 9188) hold the shaft firmly on the key surface and pull out the chuck with a jolting action or by tapping on the chuck key
3. Remove the chuck from the shaft

### 13.2 Removal of chuck on the handpiece **QUBE** long version

1. Remove tool, unscrew tip of handpiece and remove dustcap
2. Remove cap from motor and pull out the plug
3. Place the chuck key (Art. no. 4115) on the open chuck, use a number 6223 wrench to hold the motor end of the spindle
4. Remove the chuck from the shaft

After thoroughly cleaning the chuck, before replacing, it should be slightly greased on the outer surface (with grease art.-no. 51/1).  
Only fasten the chuck lightly when reinstalling! The chuck will fix itself firmly as part of its normal operation.

To optimise the precision and durability of the chuck, a tool or the factory supplied protective pin must always be fixed, even when not in use.


Never clean the handpiece with compressed air!

The required chuck keys are located in the underside of the handpiece holder.

### 13.3 Removing/changing the motor cable on the short **QUBE** handpiece

1.

2.

1. Unscrew the tip (a) of the handpiece, pull the grip casing (b), the ring (c) and cable cap (d) in the direction of the arrows

2. Detach the plug element of the motor cable from the handpiece

### 13.4. Removing/changing the motor cable on the motorhandpiece **QUBE** long version (see picture page 16)

1. Unscrew the cap (2) from the motor

2. Loosen the cable by pulling out the plug-in seal (1)

## 14. Possible error messages

The **QUBE** control unit has an intelligent control electronics system which recognises possible defects and can display them with an error code on the control unit.

If there is a defect, this is shown with a flashing red light on the controls satellites, with a flashing red SCHICK-logo on the **QUBE** control unit and with an error message with the following format "E xxyy".

Example:


Error code: E 0201 - no motor connected

Error code	Error description	Error removal
E 0201	no motor connected	Check whether handpiece is correctly connected (motor output 1 / motor output 2)
E 0801	motor overload	reduce pressure to the tool
E 0802	motor blocked	Check if chuck is closed, make sure that the shaft can be easily rotated Check with chuck key tool if chuck is properly screwed in and fastened (see also figure 13)
	cable break	replace motorcable
general error	unit flashes red, no display at the satellite	Check cabel satellite
<b>general error</b>	switching off by control satellite not possible. Bad response of Touchfields  + 	Full reset by disconnecting from mains plug. Don't touch satellite while restarting the unit.

### 14.1 Reset after error message

Error messages will disappear immediately from the display when the control level is released or the fault is eliminated. This means no waiting times for a restart.

Exception: If the table top model is used without a footswitch reset of the error messages by simultaneously press and hold the reset touch fields  +  + . A full reset can be done by pressing the touch fields  +  for 2 sec., or by disconnecting the mains supply.

If an error cannot be removed following the above description, please contact an authorised service partner or Schick directly.

## 15. Set up suggestions for the electric wax knife

The heating power of the electric wax knife can be adjusted from level 1 to 20.

The following set up suggestions are based on the Schick DesignWax and can be used as an example for individual adjustments.

Wax	Set up	Working temperature
Cervical wax lila	Level 10	approx. 125°C
Ceramic wax yellow	Level 11-12	approx. 130° - 135°C
Modelling wax grey	Level 15	approx. 140°C
Milling wax green	Level 15-17	approx. 140° - 150°C

## 16. Technical Data

	<b>QUBE</b>	<b>QUBE</b> Plus	<b>QUBE</b> Premium
Speed range:	200 - 40.000 rpm	200 - 50.000 rpm	200 - 60.000 rpm
Torque max.:	6 Ncm	7,5 Ncm	9 Ncm
Handpiece <b>QUBE</b> diameter:	max. 29 mm	Handpiece <b>QUBE</b> long version diameter:	max. 27,5 mm
length:	138 mm	length:	158 mm
weight without cable:	180g	weight without cable:	240g
weight with cable:	320g	weight with cable:	370g
cooling system:	sealed system, no ventilator	cooling system:	sealed system no ventilator
drive system:	DC-Motor, no commutator	drive system:	DC-Motor, no commutator
concentricity:	< 0,02 mm	concentricity:	< 0,02 mm
chuck:	Ø 2,35 mm	chuck:	Ø 2,35 and Ø 3,00 mm
changing of burs :	rapid change	changing of burs:	rapid change

The oscillation total value during operation is below 2,5 m/s<sup>2</sup> (EN 28662).

### Control units:

dimensions:	knee control	foot control	table top model
width:	95 mm	155 mm	155 mm
height:	203 mm	75 mm	55 mm
depth:	210 mm	220 mm	155 mm
weight:	902 g	1.820 g	555 g

Dimensions:	satellite	power supply
width:	80 mm	175 mm
height:	39 mm	35 mm
depth:	80 mm	70 mm
weight:	346 g	682 g

Operating voltage	100 - 240 V
Output	160 Watt
Overload protection	No waiting time when restarting
Safety class	III

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declare herewith, that the products

**QUBE** - consisting of  
**QUBE** - motorhandpiece 9001 or 9002 in connection with  
**QUBE** - unit 9020, 9016, 9017, 9018 and

**QUBE** Plus- consisting of  
**QUBE** Plus- motorhandpiece 9001 or 9002 in in connection with  
**QUBE** Plus- unit 9020, 9036, 9037, 9038 and

**QUBE** Premium- consisting of  
**QUBE** Premium- motorhandpiece 9001 or 9002 in connection with  
**QUBE** Premium- unit 9020, 9056, 9057, 9058

are in conformity with the following provisions of Directive:

2006/42/EG (machinery directive)  
2014/30/EU (EMV directive)  
2011/65/EU (RoHS)

Name and address of person in charge: Wolfgang Schick  
Lehenkreuzweg 12  
88433 Schemmerhofen

Schemmerhofen, March 2017

W. Schick  
Director

Subject to technical modifications

This unit complies with the current VDE (German association of electrical technicians) regulations concerning safety and suppression.  
These instructions should be readily accessible and are best kept close to the unit itself.

We would like to take this opportunity to advise you that a proper repair service and suitable qualified personnel are required for such highly developed technical equipment. SCHICK guarantees to carry out perfect repairs using original spare parts.

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Schick GmbH