


# *The Dental Think Tank* *2015 | 2016*

THINGS COULD ALWAYS BE IMPROVED.  
OR JUST BE REINVENTED.



BESSER BECHT. MADE IN OFFENBURG.

 **Becht**

ALFRED BECHT GMBH

# The little things... that made us big.

We're proud of the fact that Becht is one of the success stories that have made Germany economically strong.

Founded more than 90 years ago by man with an entrepreneurial vision, the company has been the life's work of three generations. The third owner of the business has energetically held the reins for 20 years.

It's a family business in the very best sense: small enough to retain its values while being committed to its regional roots and sustainable development.

Over the years, every new development and every new product idea has taken us forwards.

1923

Alfred Becht sets up a dental wholesale company in Offenburg under his name.

1926

The first trade mark is registered in the patent office for the Becht opaline.

1937

A success story begins: the Fräsator makes routine in the dental office clean and safe.

1950

Becht starts production of the first cleaning and disinfecting liquid: Bechtol and Carborundol.

1954

Right next door is the Black Forest. That's why from now on a stylized fir tree makes the company logo unique.

1955

In the second generation the son Hans F. Becht enters as successor to the founding father.

1959

The company grows and expands. And moves to a new building on Luisenstrasse in Offenburg.

1965

A classic with generations of dentists: Rollofix – the glass container for the hygienic storage of cotton rolls.

1973

We start the fully automated production of paper and guttapercha points for endodontic treatment.

1987

Alfred Becht becomes a private limited company.



We grow on the basis of existing business, a solid foundation always producing something new, as we increase our rates of production. And although we now do business in many countries around the world, we believe the reason we are so good is that we are based right here in southern Germany, on the edge of the Black Forest.

This commitment to our local region is reflected in our sense of responsibility to our employees, partners and customers. Reliability, fairness and social engagement are values we have always believed in, and continue to believe in. It's been that way since 1923 - and it's not about to change.



1989

The building becomes too small. The Alfred Becht GmbH moves to its new premises on Carl-Zeiss Strasse conveniently placed by the gateway to the city.

1995

On the first of January the third generation goes to work. Ulrich A. Becht takes over the company.

1995

An important milestone in hygiene, Becht introduces the Medi-Spenser system.

1998

Better safe than sorry: Becht is now certified – according to ISO 9002/40002 and directive 93/42/EC governing medical devices.

2006

The growth of the disinfection and cleaning agents have been so successful that a new state of the art manufacturing warehouse had to be built.

2008

The first batch of root canal sealer Perma Evolution is produced and put on the market.

2011

The go-ahead for the new concentrate Bechtol Premium for instrument disinfection.

2012

Room for growth: Becht expands and builds another facility for the production and storage of medical devices.

2013

The development of a new product has been completed: ProphyRepair - the first all-in-one prophylactic paste with remineralization using Hydroxylapatite.

2015

Always one step ahead. With Bechtozid sensitive we set new standards in the disinfection of surfaces.



# ***Good ideas. And the ambition to keep making them better.***



When a company still bears the name of its founder after 90 years, there must be a good reason for it. The current business owner even bears his grandfather's first name as an initial. It's the perfect

way to express what Alfred Becht GmbH stands for. We strive to achieve great things every single day through a combination of experience, skill and entrepreneurial spirit. In the way we work. And through the quality of our products.

Today, Alfred Becht GmbH is a globally active, successful medium-sized company. As innovators in dental technology we are small enough to create intelligent, innovative solutions and products that set us apart from the rest. And big enough to deliver a complete process chain: development, production, marketing and long-term service.

People in this region of Germany have a reputation for innovative and creative engineering. We combine this with hard work and the ambition to continually improve. Day in, day out, this combination of tradition, skill and determination produces smart products that make life easier for you, our customers.

We're delighted that you've chosen to place your trust in our hands. Now you can leave the good ideas to us.

Kind regards,

A handwritten signature in black ink, appearing to read 'Ulrich Becht'.

Ulrich A. Becht  
Managing Director



# Content

TREATMENT	Endodontics	4
	Prophylaxis	8
	ProfiDent Polishing Strips and Polishers	12
DENTAL OFFICE	Disinfection and Cleaning	14
	Opaline and Glass Products	23
	Cotton Wool Centers and Containers	38
	Office Administration, X-Ray Record Systems	40
	Tray Paper, Articulating Paper and Sundries	44
LABORATORY	ProfiLab Brushes and Polishers	47
	Index	54
	Contact Persons	55



Because double protection means better protection.

The perfect paste for successful endodontic treatment, **Asphaline temp** is the innovative solution for temporary fillings of gangrenous root canals. It directly assists with the healing process in the root canal itself.

**Asphaline temp** stands out because of its double effectiveness – for strong and thorough protection against infection. The paste forms a mechanical self-sterilizing barrier to prevent microorganisms from entering the root canal. Bacteria and viruses are effectively killed off at the interface between the dentin and the root filling material.

**BEST  
seller**

## Asphaline temp temporary root filling material

- ▼ encourages complete healing of infected areas
- ▼ ready-to-use paste
- ▼ radiopaque
- ▼ free of paraformaldehyde and eugenol
- ▼ adheres excellent to the dentin
- ▼ can be mixed with other water-based pastes

Tube with 7 g paste

**REF 550**



## Optimum ease of use and excellent biocompatibility.

The permanent root filling material **Perma Evolution** combines epoxy-amine technology with innovative benefits for the dental professional.

The sealer features a microbicidal effect in the liquid phase thus assuring best protection against recontamination of the root canal. The product also comes with a new reactive adhesive to make the two-component material more resilient and tight. In addition, carefully selected natural dental fillers encourage the sealer to bond with the dentin.



### Perma Evolution permanent root filling material

- ▼ free of iodoform and paraformaldehyde
- ▼ excellent biocompatibility
- ▼ microbicidal action
- ▼ outstanding dimensional stability
- ▼ strong adhesion to the dentin
- ▼ radiopaque
- ▼ easy retreatment
- ▼ suitable for all orthograde root filling techniques
- ▼ fulfils the requirements of EN ISO 6876 for dental root canal sealing materials
- ▼ working time: 15 minutes
- ▼ setting time: 24 hours

**BEST  
seller**



Automix syringe with 12 g content,  
20 mixing cannulas and 20 endo tips

**REF 555**

Refill with automix syringe with 12 g content,  
without cannulas and endo tips

**REF 555.0**

Refill with 25 mixing cannulas

**REF 555.1**

Refill with 25 endo tips

**REF 555.2**

#### Starter-Kit

Automix syringe with 3 g content,  
2 mixing cannulas and 2 endo tips

**REF 555.4**



### White Paper Points

- ISO 15-40  
Box with 200 pieces **REF 931.1**
- ISO 45-80  
Box with 120 pieces **REF 931.2**
- ISO 90-140  
Box with 120 pieces **REF 931.5**
- ISO 15-40  
Box with 500 pieces **REF 932.3**
- ISO 45-80  
Box with 300 pieces **REF 932.4**

All paper and guttapercha points are manufactured according to ISO standards and are 28 mm long. Medical Devices according to the directive 93/42/EEC.

### Color-coded Paper Points >>

- ISO 15-40  
Box with 200 pieces **REF 930.1**
- ISO 45-80  
Box with 120 pieces **REF 930.2**
- ISO 15-40  
Box with 500 pieces **REF 933.1**
- ISO 45-80  
Box with 300 pieces **REF 933.2**



**REF 930/933**



### Safety Color Paper Points

- ISO 15-40  
Box with 200 pieces **REF 937.1**
- ISO 45-80  
Box with 120 pieces **REF 937.2**



**REF 937**

### Greater Taper Paper Points >>

- GT 04**  
ISO 15-40 **REF 948.1**
- GT 06**  
ISO 15-40 **REF 949.1**
- Box with 60 pieces



All paper and guttapercha points on page 6 and 7 are available assorted or in single size packaging.





### Pink Guttapercha Points

ISO 15-40 REF 934.1

ISO 45-80 REF 934.2

ISO 90-140 REF 934.5

Box with 100 pieces



The conicity of these points correspond to the instruments for root canal treatment, e.g. for step back or crown down technique.

### Greater Taper Pink Guttapercha Points

GT 04  
ISO 15-40 REF 946.1

GT 06  
ISO 15-40 REF 947.1

Box with 60 pieces



### Color-coded Guttapercha Points

ISO 15-40 REF 935.1

ISO 45-80 REF 935.2

Box with 100 pieces



REF 935



### Safety Color Guttapercha Points

ISO 15-40 REF 936.1

ISO 45-80 REF 936.2

Box with 100 pieces



REF 936



With different pastes and types of brushes, we offer a comprehensive range of prophylactic products. Rely on professional prophylactic solutions – for all cases and all requirements. Gluten and phthalate-free.

## Prophylactic Pastes from BECHT

**Fast, easy and always reliable.**

Our prophylactic pastes are designed to prevent caries and periodontal disease. They gently remove discolorations and stubborn plaque from all tooth surfaces. Smooth and clean surfaces are reached during final polishing. Sodium fluoride promotes remineralization of the enamel and prevents inflammation of the gums. And the mint flavor makes the treatment more pleasant for your patients.



### ProfiGuard

blue, for the removal of discolorations and plaque, contains sodium fluoride (0,2%)

RDA 250

Tube with 95 g/70 ml

REF 480



### ProfiPolish

yellow, for final polishing of vital teeth and fillings, contains sodium fluoride (0,2%)

RDA 40

Tube with 95 g/74ml



REF 481

## ProphyRepair All-in-one Prophylactic Paste

Intelligent solutions for the dental office, the new ProphyRepair sets distinctive standards in prophylactic treatment making it obsolete to keep several pastes with different grit to achieve final and lasting results on tooth surfaces.

Thanks to their unique geometry the silicates in the paste are progressively activated from initial cleaning to final polishing. Additionally, ProphyRepair contains Hydroxylapatite for effective remineralization of the teeth.



### All-in-one

- ▼ cleaning
- ▼ polishing
- ▼ remineralization

**BEST**  
seller



### ProphyRepair

- ▼ mild mint flavor
- ▼ only foams up slightly
- ▼ water soluble
- ▼ gluten and phthalate-free
- ▼ contains Hydroxylapatite

white, final polishing (RDA 34)  
Tube with 95 g/75 ml

REF 486

Fresh protection for teeth and gums: our dental floss gently and effectively removes impurities and plaque from between the teeth.

The pure natural silk easily glides between the teeth without causing irritation or damaging the gums. The floss is waxed for extra protection. The mint flavor provides a pleasant fresh sensation. And to encourage remineralization of the teeth, the floss is coated with the active ingredient xylitol.



## Natural Dental Floss


in the tooth box, blistered, mint flavored with xylitol and mirror tinfoil.

waxed silk floss, 25 m

**REF 484**

waxed silk floss, 50 m

**REF 485**

 **Recommended RPM: 5.000 min-1**  
All screw-in brushes are suitable for use with the prophylaxis heads of KaVo, Sirona, W+H, Micro-Mega, etc.



## Flat Brush

with nylon filament, screw-in  
Contents: 24 pieces

**REF 206**



## Pointed Brush

with natural bristles, screw-in  
Contents: 24 pieces

**REF 207**

## Mandrels

with screw mandrel to be mounted on angled handpiece, 14 mm long



**REF 210**



with snap-on mandrel to be mounted on angled handpiece

**REF 212**



 **Recommended RPM: 12.000 min-1**



## Nylon Brush

for Paromatic-Swing  
Contents: 12 pieces/box

**REF 261.4**

Cup brushes to remove plaque and discolorations, thereby achieving an aesthetically clean tooth surface. Use ProfiGuard, ProfiPolish and ProphyRepair to reach best results. Soft quality brushes to protect the gingiva.

### Cup Brushes


4 mm, mounted on CA cups



Recommended  
RPM: max. 5.000 min-1




◀◀ with natural bristles,  
12 pieces **REF 260.21**  
100 pieces **REF 260.230**



◀◀ with nylon filament, pink  
100 pieces **REF 260.250N**



◀◀ with nylon filament, white  
12 pieces **REF 260.2N**  
100 pieces **REF 260.230N**




◀◀ with nylon filament, white  
100 pieces **REF 260.260N**

### Prophy Cups latex-free

Prophy cups with lamellae for a thorough cleaning and polishing. The rubber cups are available in eye catching colors and different shapes and hardness. The underside is formed to prevent paste from penetrating the prophylaxis head. The different mounting of the polishers offers almost unlimited versatility.



**BLUE:** soft, for the removal of plaque and for polishing work, recommended for use near the sulcus for protection of the soft tissue.  
**YELLOW:** hard, for removal of heavy plaque.



Color	Mounted on CA	Screw-in cups	Snap-on cups
BLUE	REF 204.1	REF 200.1	REF 202.1
YELLOW	REF 204.2	REF 200.2	REF 202.2

Contents: 24 pieces

Recommended  
RPM: max. 10.000 min-1

Becht brushes and rubber cups for prophylaxis are medical devices, class IIa. For single use only. All packing units are supplied unsterile. The low friction coefficient of the material prevents the enamel to heat up.

# ProfiDent Polishing Strips and Polishers



## Amalgam Polishing Strips

Highly flexible strips for polishing the interproximal area of amalgam fillings. The transparent base with its blue carborundum particles makes for consistently good visibility during restoration. Efficient removal of excess material with no risk of injury to the gums. With easy-grip ends and uncoated window for easy insertion between the teeth. Coated on one side.

Grit: 80 µ. Length: 150 mm

4 mm wide

6 mm wide

Contents: 12 pieces

REF 9045.04

REF 9045.06

## Steel Separating Strips

Flexible strips for removal of excess material of amalgam fillings. With galvanically bonded aluminium oxide. Coated on one side. Suitable for separating teeth. The strips can be cut to the required length.

Grit: 80 µ, Length: 135 mm

4 mm wide

6 mm wide

Contents: 12 pieces

REF 9042.04

REF 9042.06



## Polishing/High-shine Polishing of Precious Alloys/Amalgam/Composites



REF	242.1	242.2	242.3
Application	Polishing of precious alloys/amalgam/composites	Polishing of precious alloys/amalgam/composites	Polishing of precious alloys/amalgam/composites
RPM	7.000 - 12.000 min-1	7.000 - 12.000 min-1	7.000 - 12.000 min-1
Product description	Silicone polisher, pointed, medium	Silicone polisher, small point, medium	Silicone polisher, cup, medium
Size	4 x 12 mm	3 x 6 mm	6,5 x 9,5 mm




REF	243.1	243.2	243.3
Application	High-shine polishing of precious alloys/amalgam/composites	High-shine polishing of precious alloys/amalgam/composites	High-shine polishing of precious alloys/amalgam/composites
RPM	7.000 - 12.000 min-1	7.000 - 12.000 min-1	7.000 - 12.000 min-1
Product description	Silicone polisher, pointed, fine	Silicone polisher, small point, fine	Silicone polisher, cup, fine
Size	4 x 12 mm	3 x 6 mm	6,5 x 9,5 mm

### Starter Kit silicone polishers


with one polisher each 242.1 – 243.3 (6 pieces)

REF 243.SET

## Pre-polishing/Polishing/High-shine Polishing of Ceramics/Composites



REF	240.1	240.6	240.7	240.8
Application	Pre-polishing of ceramics/composites	Pre-polishing of ceramics/composites	Pre-polishing of ceramics/composites	Pre-polishing of ceramics/composites
RPM	7.000 - 12.000 min <sup>-1</sup>	7.000 - 12.000 min <sup>-1</sup>	7.000 - 12.000 min <sup>-1</sup>	7.000 - 12.000 min <sup>-1</sup>
Product description	Silicone polisher, flame, medium	Silicone polisher, lens, medium	Silicone polisher, pointed, medium	Silicone polisher, cup, medium
Size	5,5 x 15 mm	11 x 2,5 mm	5 x 10 mm	7 x 10 mm



REF	226.1	226.2	227.1	227.2
Application	Polishing of ceramics	Polishing of ceramics	High-shine polishing of ceramics	High-shine polishing of ceramics
RPM	7.000 - 15.000 min <sup>-1</sup>	7.000 - 15.000 min <sup>-1</sup>	7.000 - 15.000 min <sup>-1</sup>	7.000 - 15.000 min <sup>-1</sup>
Product description	Silicone polisher, pointed, medium	Silicone polisher, cup, medium	Silicone polisher, pointed, fine	Silicone polisher, cup, fine
Size	5 x 10 mm	7 x 10 mm	5 x 10 mm	7 x 10 mm

### Starter Kit silicone polishers

with one polisher each 240.1 – 240.8 (4 pieces)

**REF 240.SET**

with one polisher each 240.7 – 227.2 (6 pieces)

**REF 227.SET**

## Pre-polishing/Polishing of Composites



REF	228.1	228.2	228.3
Application	Pre-polishing of composites	Pre-polishing of composites	Pre-polishing of composites
RPM	7.000 - 12.000 min <sup>-1</sup>	7.000 - 12.000 min <sup>-1</sup>	7.000 - 12.000 min <sup>-1</sup>
Product description	Polisher of Flex-material, pointed, medium	Polisher of Flex-material, cup, medium	Polisher of Flex-material, lens, medium
Size	5 x 10 mm	6,5 x 9,5 mm	10 x 2,5 mm



REF	229.1	229.2	229.3
Application	High-shine polishing of composites	High-shine polishing of composites	High-shine polishing of composites
RPM	5.000 - 12.000 min <sup>-1</sup>	5.000 - 12.000 min <sup>-1</sup>	5.000 - 12.000 min <sup>-1</sup>
Product description	Polisher of Flex-material, pointed, fine	Polisher of Flex-material, cup, fine	Polisher of Flex-material, lens, fine
Size	5 x 10 mm	6,5 x 9,5 mm	10 x 2,5 mm

### Starter Kit silicone polishers

with one polisher each 228.1 – 229.3 (6 pieces)

**REF 229.SET**



**ProfiDent**

polishers can easily be reprocessed in the autoclave (up to 134 °C). Sterilization must not be performed in a chemical autoclave. For disinfection we recommend our tried and tested instrument disinfectants **Bechtol Premium/Bechtol plus**. Observe instructions of use. The ProfiDent line of polishers are medical devices according to the directive 93/42/EEC. All polishers are supplied for **contra-angle hand-piece, unsterile, in boxes of 6 pieces.**



The practical everyday range for regular hygiene procedures in dental and podiatry offices. Products that users are happy to recommend.



### Bechtol Premium Concentrate

Next generation disinfecting and cleaning solution for surgical instruments

- ▼ highly effective against bacteria (including MRSA) and fungi in 15 minutes in a 2% solution. Guarantees active protection against *M. terrae* in 60 minutes in a 1% solution
- ▼ tested in accordance with the directives of the DGHM (Deutsche Gesellschaft für Hygiene und Mikrobiologie) of 01.09.2001 under high organic load

- ▼ inactivates enveloped viruses (HBV/HCV/HIV/BVDV/Vaccinia) and non-enveloped viruses (Adeno/Polio/HPV/Rota and Noro) in a 3% solution in 15 minutes (EN 14476)
- ▼ contains no phenols or aldehydes
- ▼ suitable for use in an ultrasound device in 5 minutes in a 3% solution
- ▼ excellent compatibility of materials

1-Liter dispenser bottle, concentrate, reusable

**REF 536**

5-Liter canister, concentrate

**REF 536.5**

**BEST**  
seller





## Bechtol CLASSIC Concentrate

For the disinfection and cleaning of surgical instruments. Gentle on materials. High-performance. Economical.

- ▼ effective against bacteria (incl. MRSA) and fungi in 15 minutes in a 3% solution (EN 13624/13727)
- ▼ guarantees active protection against tuberculosis (*M. terrae*) in 30 minutes in a 3% solution under high organic load (EN 14348)
- ▼ inactivates enveloped viruses (HBV/HCV/HIV/BVDV/Vaccinia) and non-enveloped viruses (SV40/Noro) in 15 minutes in a 3% solution
- ▼ contains no phenols or aldehydes
- ▼ biodegradable
- ▼ suitable for use in an ultrasound device in 5 minutes in a 4% solution
- ▼ excellent compatibility on materials


1-Liter dispenser bottle, concentrate, reusable

**REF 534**

5-Liter canister, concentrate

**REF 534.5**



 Rotating instruments can hygienically and economically be cleaned in our bur disinfectors "Fräsator" REF 630 and REF 631 (description on page 30).



## Bechtol plus Ready-to-use Solution

The highly effective disinfecting and cleaning solution for rotating instruments

- ▼ effective against bacteria (incl. MRSA) and fungi in 5 minutes
- ▼ kills spores in 60 minutes
- ▼ tested in accordance with the directives of the DGHM (Deutsche Gesellschaft für Hygiene und Mikrobiologie) of 01.09.2001

- ▼ inactivates HBV/HCV/HIV/HPV/Adeno/Polio viruses in 5 minutes (EN 14476)
- ▼ contains no phenols or aldehydes
- ▼ with fresh citrus scent
- ▼ biodegradable

1-Liter bottle ready-to-use solution

**REF 539**

5-Liter canister ready-to-use solution

**REF 539.5**



### Bechtol Ready-to-use Solution

The non-corrosive disinfecting and cleaning bath for burs, diamonds and root canal instruments

- ▼ kills bacteria and fungi in 30 minutes and inactivates HBV/HCV/HIV viruses in 60 minutes
- ▼ tested in accordance with the directives DIN EN 13727 and DIN EN 13624
- ▼ disinfection, cleaning and lubricating in one step
- ▼ gentle on materials even during extended soaking, e.g. over the weekend
- ▼ non-corrosive
- ▼ contains no phenols or aldehydes
- ▼ biodegradable
- ▼ product known for its reliability
- ▼ suitable for use in an ultrasound device

1-Liter bottle ready-to-use solution

**REF 543**

5-Liter canister ready-to-use solution

**REF 543.5**

### Becht Carborundol Ready-to-use Solution



For gentle disinfection and cleaning of trimmers, polishing stones and tips

- ▼ effective against bacteria and yeast in 60 minutes, inactivates viruses in minimal time
- ▼ tested in accordance with the directives DIN EN 1040 and DIN EN 1275
- ▼ cleans gently without harming tip bonding
- ▼ non-corrosive
- ▼ contains no phenols or aldehydes
- ▼ biodegradable

1-Liter bottle ready-to-use solution

**REF 562**

5-Liter canister ready-to-use solution

**REF 562.5**





**!** Ensure maximum safety for yourself and your patients by using the correct dilution. You will find a dilution table on our website [www.becht-online.de](http://www.becht-online.de)

## Bechtid plus Concentrate

The highly effective disinfecting and cleaning solution for use on large surfaces, e.g. floors

- ▼ effective against bacteria (incl. MRSA) and fungi in 15 minutes in a 2% solution (EN 1040/1276)
- ▼ tested in accordance with the directives of the DGHM (Deutsche Gesellschaft für Hygiene und Mikrobiologie) of 01.09.2001
- ▼ inactivates HBV/HCV/HIV/Vaccinia/HPV/Adeno/Polio viruses in 15 minutes in a 5% solution

- ▼ disinfects, cleans and protects in just one application on all wet proved surfaces
- ▼ contains no phenols or aldehydes
- ▼ biodegradable

1-Liter dispenser bottle, concentrate, reusable

**REF 546**

5-Liter canister, concentrate

**REF 546.5**



## Bechtozid soft Ready-to-use Solution

Foam-spray for the disinfection and cleaning of all alcohol sensitive medical devices, natural and synthetic leather

- ▼ highly effective against bacteria (incl. MRSA), fungi and viruses (HBV/HCV/HIV/BVDV/Vaccinia) in 1 minute
- ▼ tested in accordance with the directives of the DGHM (Deutsche Gesellschaft für Hygiene und Mikrobiologie) of 01.09.2001
- ▼ alcohol-free and without aerosol
- ▼ free of perfumes and coloring agents for best compatibility
- ▼ duo spray head (foam or spray)

1-Liter bottle, ready-to-use solution, reusable, duo spray head

**REF 537**

5-Liter canister, ready-to-use solution

**REF 537.5**

10-Liter canister, ready-to-use solution

**REF 537.10**

Duo spray head for 1-Liter bottle

**REF 537.SP**



## Bechtozid Premium Ready-to-use Solution



For the disinfection and cleaning of high-risk surfaces of medical devices, e.g. patient treatment chairs, unit handles and handpieces

- ▼ highly effective against bacteria (incl. MRSA) and fungi in 30 sec
- ▼ tuberculocidal (M. terrae) in 30 seconds (EN 14348)
- ▼ inactivates HBV/HCV/HIV/HPV/BVDV/Vaccinia viruses in 30 sec (EN 14476)
- ▼ inactivates Adeno in 1 minute and Polio in 2 minutes (EN 14476)
- ▼ effective against H1N1-viruses (swine flu) in 30 sec (EN 14476)
- ▼ tested in accordance with the directives of the DGHM (Deutsche Gesellschaft für Hygiene und Mikrobiologie) of 01.09.2001 under high organic load
- ▼ contains no phenols or aldehydes
- ▼ fast residue-free drying
- ▼ gentle to all alcohol resistant surfaces
- ▼ non-corrosive
- ▼ perfume-free or with fresh citrus scent



1-Liter bottle, ready-to-use solution, perfume-free, reusable, without spray head

**REF 561PF**

1-Liter bottle, ready-to-use solution, with citrus scent, reusable, without spray head

**REF 561**

Spray head for 1-Liter bottle

**REF 544.SP**



5-Liter canister, ready-to-use solution, perfume-free

**REF 561.5PF**

5-Liter canister, ready-to-use solution, with citrus scent

**REF 561.5**

10-Liter canister, ready-to-use solution, perfume-free

**REF 561.10PF**

10-Liter canister, ready-to-use solution, with citrus scent

**REF 561.10**



## Bechtozid plus Ready-to-use Solution

The rapid surface disinfectant for optimal patient flow

- ▼ effective against bacteria (incl. MRSA) and fungi in 1 minute
- ▼ tested in accordance with the directives of the DGHM (Deutsche Gesellschaft für Hygiene und Mikrobiologie) of 01.09.2001 under high organic load
- ▼ inactivates HBV/HCV/HIV/BVDV/Vaccinia viruses in just 1 minute
- ▼ contains no phenols or aldehydes
- ▼ gentle to all alcohol resistant surfaces (Röhm test 12/2006)
- ▼ neutral or with fresh citrus scent
- ▼ biodegradable

1-Liter bottle, ready-to-use solution, neutral, reusable, without spray head

**REF 544N**

1-Liter bottle, ready-to-use solution, with citrus scent, reusable, without spray head

**REF 544**

Spray head for 1-Liter bottle

**REF 544.SP**

5-Liter canister, ready-to-use solution, with citrus scent

**REF 544.5**

10-Liter canister, ready-to-use solution, neutral

**REF 544.10N**

10-Liter canister, ready-to-use solution, with citrus scent

**REF 544.10**



## Bechtozid sensitive Ready-to-use Solution

QAC-free surface disinfecting and cleaning spray for medical devices, e.g. trays, hand pieces.  
High-performance. Gentle on materials. Economical.

- ▼ effective against bacteria (incl. MRSA) and fungi in 1 minute (EN 13624/13727)
- ▼ guarantees active protection against tuberculosis (M. terrae) in 5 minutes under high organic load (EN 14348)
- ▼ inactivates HBV/HCV/HIV/BVDV/Vaccinia-viruses in 1 minute
- ▼ avoids irritation of the respiratory system
- ▼ low contents of alcohol (< 30%)

- ▼ for use on acrylic surfaces
- ▼ neutral or with fresh lemon scent

- |   |                    |
|---|--------------------|
| 1-Liter bottle, ready-to-use solution, lemon scent      | <b>REF 563L</b>    |
| 1-Liter bottle, ready-to-use solution, neutral scent    | <b>REF 563N</b>    |
| Spray head for 1-Liter bottle                           | <b>REF 544.SP</b>  |
| 5-Liter canister, ready-to-use solution, lemon scent    | <b>REF 563.5L</b>  |
| 5-Liter canister, ready-to-use solution, neutral scent  | <b>REF 563.5N</b>  |
| 10-Liter canister, ready-to-use solution, lemon scent   | <b>REF 563.10L</b> |
| 10-Liter canister, ready-to-use solution, neutral scent | <b>REF 563.10N</b> |



## Bechtifix plus

For disinfecting and cleaning of medical devices, e.g. treatment chairs, trays and handles in dental, medical and chiropody offices

- ▼ effective against bacteria (incl. MRSA) and fungi in 1 minute under high organic load (EN 13624/13727).
- ▼ **successfully tested in the 4-field test acc. to EN 16615**
- ▼ inactivates HBV/HCV/HIV/BVDV/Vaccinia-viruses in 1 minute (EN 14476)
- ▼ **new quality of tissues:** particularly soft, absorbent and tear-resistant
- ▼ high volume of release of disinfecting solution on surface
- ▼ gentle to all alcohol resistant surfaces
- ▼ suitable for use on perspex (Röhm-test)
- ▼ with fresh lemon scent
- ▼ size of tissue: 14 x 20 cm

Dispenser box with 100 tissues

**REF 548.100L**

Refill with 100 tissues

**REF 548.100LNF**



Available in summer 2015



Available in summer 2015

## Bechtifix alcohol-free

Tissues for the disinfection and cleaning of all alcohol sensitive medical devices, natural and synthetic leather



- ▼ effective against bacteria (incl. MRSA) and fungi in 1 minute under high organic load (EN 13624/13727)
- ▼ **successfully tested in the 4-field test acc. to EN 16615**
- ▼ inactivates HBV/HCV/HIV/BVDV/Vaccinia-viruses in 1 minute (EN 14476)
- ▼ particularly soft, absorbent and tear-resistant non-woven fabric
- ▼ high volume of release of disinfecting solution on surface
- ▼ free of perfumes and coloring agents for best compatibility
- ▼ size of tissue: 14 x 20 cm

Dispenser box with 100 tissues

**REF 549.100L**

Refill with 100 tissues

**REF 549.100NF**

## Asphalin B Tablets

To keep sterilized dental and surgical instruments sterile until reuse

- ▼ highly efficient and most economic in use (14 days approx. / tablet)
- ▼ specially recommended for root canal instruments stored in Asphalin trays
- ▼ easy to handle cell pack packing

Cell pack packing with 60 tablets

REF 554

Cell pack packing with 10 tablets

REF 554.1



## AirFresh med

Air cleaner



- ▼ cleans the air from dust and germs
- ▼ specially recommended for use in dental clinics, hospitals and nursing homes
- ▼ odour-absorbing
- ▼ lasting fresh scent
- ▼ economical use
- ▼ contains propane/butane as propellant
- ▼ biodegradable acc. to OECD guidelines

75 ml spray can

REF 590

## AquaHy Concentrate

Stabilizer liquid for spray devices used in podiatry

- ▼ reduces germs within a few minutes
- ▼ without scent
- ▼ easy to use
- ▼ economical (concentration 1:50)

Bottle with 100 ml concentrate

REF 560






## Bur and Diamond Cleansing Brushes

with plastic handle, can be disinfected

with brass wire, thickness 0,08 mm  
Spare brass wire

**REF 145**  
**REF 145.1**



with stainless steel wire, thickness 0,08 mm  
Spare stainless steel wire

**REF 147**  
**REF 147.1**



## Bur Brush

adjustable, with brass wire, thickness, 0,10 mm  
with stainless adjustable belt

**REF 144**



## Disinfection tub

capacity 3 liter approx., with drain inset and lid,  
of plastic material

Size: 315 x 206 x 125 mm

**REF 566**

## Sundries



Dispenser pump for 1-Liter bottle  
Dispenser pump for 5 Liter canister  
Dispenser pump for 10 Liter canister

**REF 567.1**  
**REF 567.5**  
**REF 567.10**



Pouring spot for 5 and 10-Liter canister  
Canister key for 5 and 10-Liter canister

**REF 567.51**  
**REF 567.61**







Dispenser bottles permit the economical, hygienic and exact dosage of medicaments. The Medi-Spensers come in different colors to avoid confusion of content. The autoclavable 10 ml bottles are supplied complete with membranes, combi stand and dappen dishes.



### Medi-Spenser Combination System

complete with Medi-Spensers, silicone membranes, dappen dishes, combi stand

Set of 5, available in the basic colors crystal, yellow, rose, green, blue

**REF 796.0**

Set of 3, available in the basic colors crystal, yellow, rose

**REF 796.1**



### Combi Stands

of chrome-nickel steel, stainless

for 5 Medi-Spensers, 5 dappen dishes  
Size: 198 x 98 x 33 mm

**REF 395**

for 3 Medi-Spensers, 3 dappen dishes  
Size: 122 x 98 x 33 mm

**REF 395.1**

### Medi-Spensers

with silicone membranes

Set of 5, available in the colors crystal, yellow, rose, green, blue

**REF 795.0**

Set of 3, available in the colors crystal, yellow, rose

**REF 795.1**

### Medi-Spensers

Individual colors, with silicone membranes

crystal  
yellow  
rose  
green  
blue  
brown

**REF 795.3**

**REF 795.4**

**REF 795.5**

**REF 795.6**

**REF 795.7**

**REF 795.8**

### Silicone Membranes

for Medi-Spensers, autoclavable

Set of 5 membranes

**REF 795.10**

## Dappen Dishes



of opaline, autoclavable. They come with an easy-grip, corrugated outer surface and can be used on both ends for different purposes, e.g. for the preparation of cement and filling materials as well as for the economical and easy loading of the amalgam carrier. Approx. 2-3 ml capacity.

Size: 32 mm in dia., 30 mm high

REF 689



## Dappen Dishes

of colored glass, decagonal, with easy-grip surface. Both ends with different shapes. They are suitable for the preparation of cement and filling materials. The dappen dishes are also supplied individually.

Not autoclavable. Approx. 2-3 ml capacity.

Size: 32 mm in dia., 30 mm high.

crystal, yellow, rose, green, blue  
Set of 5 glasses, assorted

REF 783.7

crystal, yellow, rose  
Set of 3 glasses, assorted

REF 783.5

Individual colors

crystal  
yellow  
rose  
green  
blue  
brown

REF 783  
REF 784  
REF 785  
REF 786  
REF 787  
REF 788

For information on proper reprocessing see  
[www.becht-online.de](http://www.becht-online.de)



## Medicament Bottles



of colored glass

crystal, yellow, rose, green, blue

Set of 5 glasses, assorted

REF 743.7

	Bottle with lid
crystal	REF 743
yellow	REF 744
rose	REF 745
green	REF 746
blue	REF 747
brown	REF 748

Lids
REF 743.2
REF 744.2
REF 745.2
REF 746.2
REF 747.2
REF 748.2





## Medicament Bottles

of high quality opaline

autoclavable, with color rings on the shoulder of the base and lower brim of the cover.

Set of 5 colors, assorted

**REF 699.0**

Color ring	Bottle with lid
yellow	<b>REF 699.02</b>
red	<b>REF 699.03</b>
green	<b>REF 699.04</b>
blue	<b>REF 699.05</b>
black	<b>REF 699.06</b>

Lid
<b>REF 699.22</b>
<b>REF 699.23</b>
<b>REF 699.24</b>
<b>REF 699.25</b>
<b>REF 699.26</b>

## Medicament Bottles

of high quality opaline with special color coding system for storage in drawers

autoclavable, with color rings on the upper brim of the base and lower brim of the cover.

Set of 5 colors, assorted

**REF 699.1**

Color ring	Bottle with lid
yellow	<b>REF 699.12</b>
red	<b>REF 699.13</b>
green	<b>REF 699.14</b>
blue	<b>REF 699.15</b>
black	<b>REF 699.16</b>

Lid
<b>REF 699.22</b>
<b>REF 699.23</b>
<b>REF 699.24</b>
<b>REF 699.25</b>
<b>REF 699.26</b>



## Medicament Bottles

of Simplex opaline

autoclavable, with color rings on the shoulder of the base and lower brim of the cover.

Set of 5 colors, assorted

**REF 698.0**

Color ring	Bottle with lid
yellow	<b>REF 698.02</b>
red	<b>REF 698.03</b>
green	<b>REF 698.04</b>
blue	<b>REF 698.05</b>
black	<b>REF 698.06</b>

Lid
<b>REF 698.22</b>
<b>REF 698.23</b>
<b>REF 698.24</b>
<b>REF 698.25</b>
<b>REF 698.26</b>

Bottles for medicaments. The coding system with rings in different colors ensures safe identification of its content. Lid with easy-grip cover and airtight closure. The well rounded interior permits the contents to be used up to the last drop. Size: 32 mm in dia. at base, 42 mm high, approx. 3 ml capacity.



### Bottle Stand

of opaline, autoclavable. The curved bottle stand holds 5 dappen dishes or covered bottles.

Size: 215 x 24 mm

REF 623

### Bottle Stand

of opaline, autoclavable. The straight bottle stand holds 5 dappen dishes or covered bottles.

Size: 207 x 47 x 20 mm



REF 600



### Bottle Stand

of chrome-nickel steel, stainless, matt surface, straight.  
for 5 dappen dishes or covered bottles

Size: 198 x 44 x 17 mm

REF 601

for 3 dappen dishes or covered bottles

Size: 121 x 44 x 17 mm

REF 601.1



### Paper Tissues

of crepe paper, highly absorbent, to collect liquids, in use with dappen dishes or covered bottles stored in bottle stands or dental drawers.

Size: 213 x 75 mm

for 5 dappen dishes or covered bottles

Pack of 250 pieces



REF 33



## Standard Trays and Containers

Autoclavable. Our organization aids are available in different sizes and with compartments for the orderly placing of instruments and material for reprocessing and storage. Choose according to your individual needs.



**!** To comply with official rules always use the trays and containers with lid (see page 29) to protect your instruments against impurities and contamination. Only use disinfecting liquids recommended by the manufacturer.

### 1 Standard Tray

of opaline, all-purpose tray for the hygienic placing of dental forceps, elevators, needle holders, scissors, mirrors and tweezers.

Size: 200 x 100 x 22 mm

REF 1

### 2 Standard Tray

of opaline, all-purpose tray for the hygienic placing of 8 single-ended instruments.

Size: 200 x 100 x 30 mm

REF 2

### 3 Standard Tray

of opaline, for the orderly placing of small instruments in 6 compartments.

Size: 200 x 100 x 22 mm

REF 3

### 4 Standard Tray

of opaline, glass knobs will keep up to 9 double-ended instruments in their place.

Size: 200 x 100 x 22 mm

REF 4

### 5 Standard Tray

of opaline, for hand instruments and separating strips. Recesses in the partitions make it easy to pick up the instruments from the 3 compartments.

Size: 200 x 100 x 22 mm

REF 6

### 6 Standard Tray

of opaline, instruments of various length can be arranged in 6 compartments.

Size: 200 x 100 x 22 mm

REF 8

### 7 Standard Tray

of opaline, to rest or exchange burs. This tray makes use of even small spaces in drawers.

Size: 100 x 50 x 22 mm

1/4 standard size

REF 14

### 8 Standard Tray

of opaline, for burs and cotton pellets. This tray makes use of even small spaces in drawers.

Size: 100 x 100 x 22 mm

1/2 standard size

REF 15

### 9 Antifekt Instrument Rest

of opaline, for single-ended instruments.

Size: 120 x 30 x 25 mm

REF 627

### 10 Porcelain Tray

autoclavable, for use with small instruments, designed for multi-purpose

Size: 70 x 40 x 18 mm

REF 1718



### Crystal Glass Instrument Tray

Size: 120 x 60 x 40 mm

REF 760

Size: 200 x 100 x 50 mm, autoclavable

REF 761

Size: 200 x 100 x 50 mm

REF 762

### Opaline Instrument Tray

Size: 200 x 100 x 55 mm, autoclavable

REF 656

### Denture Bowl

of opaline, autoclavable, for universal use

Size: 110 mm in dia., 35 mm high



REF 668



### UNIKAP Closure

of opaline, autoclavable to cover the outlet of spittoon bowls, suitable for all brands

Size: 65 mm in dia., 41 mm high

REF 639

Size: 60 mm in dia., 41 mm high

REF 639.1



## 1 Stainless Steel Lid

with knob, but without rubber packing and counter glass pane

- REF 360.2                      Tray 760
- REF 362.2                      Tray 761
- REF 362.2                      Tray 762
- REF 362.2                      Standard Tray REF 1-8

## 2 Stainless Steel Lid


without knob, rubber packing and counter glass pane

- REF 362.3                      Tray 761
- REF 362.3                      Tray 762
- REF 362.3                      Standard Tray REF 1-8



### Lids for Trays

of chrome-nickel steel, with matt surface. Choose the rubber packing and counter glass pane type to assure airtight closure.



### Stainless Steel Lid

with knob, rubber packing and counter glass pane

- REF 362.1                      Tray 761
- REF 362.1                      Tray 762
- REF 362.1                      Standard Tray REF 1



### Waste Receiver

of porcelain, autoclavable

Size: 72 mm in dia., 115 mm high

REF 1785

## Disinfecting Jar

of opaline, autoclavable. It has a broad base for better stability and a reinforced brim to avoid chipping when hit by instruments.

Size: 92 mm in dia., 122 mm high

REF 648



## Bur Disinfector "Fräsator"


The burs are placed into the disinfection liquid through the slit in the cover. The sieve-insert can be lifted to allow the liquid to drain. The "Fräsator" has proven its usefulness for years and has become an indispensable asset of state-of-the-art disinfection methods. The bur disinfector "Fräsator" is available in two models.



**BEST  
seller**



REF 630.1/631.1



REF 630.2/631.2



REF 630.3/631.3

### 1 Fräsator

of opaline, the "Fräsator" is autoclavable and comes complete with slit cover and sieve-insert for all bur types and root canal instruments.

Size: 96/80 mm in dia., 65 mm high

**REF 630**

#### Spare parts:

Base  
Sieve-insert (plastic)  
Cover with slit  
Closure rod (without illustration)

**REF 630.1**  
**REF 630.2**  
**REF 630.3**  
**REF 630.4**

### 2 Mini-Fräsator

of opaline, the "Mini-Fräsator" is autoclavable and comes complete with slit cover and sieve-insert. This model is particularly suitable for turbine burs and angled instruments.

Size: 78/60 mm in dia., 66 mm high

**REF 631**

#### Spare parts:

Base  
Sieve-insert (plastic)  
Cover with slit  
Closure rod (without illustration)

**REF 631.1**  
**REF 631.2**  
**REF 631.3**  
**REF 631.4**

### Bur Tweezers

overlapping claws of the tweezers make it easy to safely pick up and lift the instruments from the disinfection liquid without the risk of bending or breaking root canal instruments.

**REF 910**





### 1 Waste Receiver

free standing model. The opaline waste receiver is autoclavable and complete with exchangeable lid no. 653 with star-shaped opening. Used cotton pads and pellets can easily be removed.

Size: 75 x 72 mm

REF 652

### 2 Waste Receiver

of opaline. The waste receiver is autoclavable and complete with exchangeable lid no. 653 with star-shaped opening.

Size: 80 x 66 mm

REF 663

### 1 Cotton Wool Container



free standing model. The opaline container is autoclavable and complete with exchangeable lid no. 655 and coil spring no. 327 for pushing up the cotton wool.

Size: 75 x 72 mm

REF 654

### 2 Cotton Wool Container

of opaline. The container is autoclavable and complete with exchangeable lid no. 655 with circular opening and with coil spring no. 327 for pushing up the cotton wool.

Size: 80 x 66 mm

REF 664



### 1 Pellet Container

free standing model. The opaline container is autoclavable and complete with metal insert and grid no. 326 and coil spring no. 327 for pushing up the pellets.

Size: 75 x 72 mm

REF 651

### 2 Pellet Container

of opaline. The container is autoclavable and complete with metal insert and grid no. 326 and coil spring no. 327 for pushing up the pellets.

Size: 80 x 66 mm

REF 661

### Glass Jars for Cotton Wool, Pellets and Dabs

the contents of the jars can be kept sterile with ground in or ground on glass covers.

with ground in glass cover, autoclavable

Size: 65 mm in dia., 40 mm high

REF 727

with ground on glass cover, not autoclavable

Size: 77/92 mm in dia., 79 mm high

REF 728



### Glass Jars for Cotton Wool

autoclavable, with overlapping glass cover.

Size: 100 mm in dia., 50 mm high

REF 734

Size: 100 mm in dia., 100 mm high

REF 735

Size: 120 mm in dia., 120 mm high

REF 736

### Glass Jars for Cotton Wool

autoclavable, with ground in glass cover to keep cotton wool sterile.

Size: 100 mm in dia., 100 mm high

Size: 120 mm in dia., 120 mm high

REF 735.1

REF 736.1



For information on proper reprocessing see [www.becht-online.de](http://www.becht-online.de)



## Asphalin Trays

the different sizes of the compartments allow the individual filling with root canal instruments, paper points and guttapercha points. Airtight closure.

Size: 135 mm in dia., 45 mm high

of opaline, 8 compartments, autoclavable

**REF 633**

of crystal glass, 8 compartments, not autoclavable

**REF 700**

## Bur Tweezers

overlapping claws ensure the safe and easy picking up of all small objects thus preventing the bending or breaking of the root canal instruments.

**REF 910**

## Petri Dishes

of crystal glass, autoclavable and resistant to sudden changes of temperature. Available with different partitions and in several sizes for the orderly keeping of small parts.



## Petri Dish

without partitions

Size: 60 mm in dia., approx. 23 mm high

**REF 791**

## Petri Dish

without partitions

Size: 80 mm in dia., approx. 23 mm high



**REF 792**





### Petri Dish

without partitions  
Size: 100 mm in dia.,  
approx. 20 mm high


REF 703

### Petri Dish



with two compartments  
Size: 100 mm in dia.,  
approx. 23 mm high

REF 704



### Petri Dish

with three compartments  
Size: 100 mm in dia.,  
approx. 23 mm high


REF 705

### Petri Dish



with four compartments  
Size: 100 mm in dia.,  
approx. 23 mm high

REF 706



**BEST**  
seller



### Alcohol Lamp

of clear glass, complete with ground on cover, burner and wick.  
Capacity: approx. 75 ml.

Size: 64 mm in dia., 90 mm high

REF 742

#### Spare parts:

Cap for alcohol lamp

REF 742.1

Plaited wick, 6 mm in dia.,  
by the metre

REF 458

Burner tube with cork

REF 468



### Alcohol Lamp

of opaline, complete with ground on cover, burner and wick.  
Capacity: approx. 75 ml.

Size: 64 mm in dia., 90 mm high

REF 640

#### Spare parts:

Cap for alcohol lamp

REF 640.1

Plaited wick, 6 mm in dia.,  
by the metre

REF 458

Burner tube with cork

REF 468



### Alcohol Lamp

of clear glass, complete with screw cap and adjustable wick.  
Capacity: approx. 100 ml.

Size: 80 mm in dia., 82 mm high

REF 755

#### Spare parts:

Plaited wick, 5 mm in dia., by the metre

REF 457

Screw cap

REF 469



### Alcohol Lamp

of opaline, complete with screw cap and adjustable wick.  
Capacity: approx. 100 ml.

Size: 65 mm in dia., 92 mm high

REF 643

#### Spare parts:

Plaited wick, 5 mm in dia.,  
by the metre

REF 457

Screw cap

REF 469





## Alcohol Lamp

of clear glass, with lateral inlet, complete with ground on cap, burner and wick. Capacity: approx. 90 ml.

Size: 79 mm in dia., 123 mm high

REF 741

### Spare parts:

Plaited wick, 6 mm in dia., by the metre

REF 458

Burner tube with cork

REF 468

## Alcohol Lamp



of stainless steel, complete with screw cap and adjustable wick. Capacity: approx. 60 ml.

Size: 70 mm in dia., 64 mm high

REF 329

### Spare parts:

Flat wick, 12 mm, by the metre

REF 460



## Reusable rinsing jars

are not only practical for your patients, they also reduce the waste volume in your dental office, thus contributing to the protection of the environment.



## Rinsing Jars "ÖKOPAL"

of opaline, stackable, autoclavable. The rinsing jars are dishwasher-safe and hardened to withstand easy breaking. Space-saving!

**BEST  
seller**

Size: upper end: 72 mm in dia., base: 45 mm in dia., 70 mm high.

Weight: 116 g

Pack of 12 jars

REF 650

## Rinsing Jars



cylindrical jars of crystal glass, with melted brim, autoclavable and dishwasher-safe. Size: 64 mm in dia., 86 mm high.

Weight: approx. 65 g

Pack of 12 jars

REF 658



**Glass Mixing Blocks and Slabs**

Practical for mixing cements. Ground or polished surface on all sides. Bevelled edges. Large mixing area. The thickness of the glass prevents an increase of temperature and thus the premature setting of the cement.



**1 Glass Mixing Block**

one side ground to avoid sliding on the desk surface.

Size: 150 x 75 mm, approx. 18-20 mm thick

REF 722

**2 Glass Mixing Block**

both sides polished.

Size: 150 x 75 mm, approx. 18-20 mm thick

REF 723

**1 Glass Slabs**



one side ground to avoid sliding on the desk surface.

Size: 95 x 70 mm, approx. 5-6 mm thick


REF 724

**2 Glass Slabs**

both sides polished.

Size: 95 x 70 mm, approx. 5-6 mm thick

REF 725



The glass slabs are supplied in packs of 12 each.

## Cotton Wool Center

complete with dividers for universal use of cotton pellets and pads.

smoked  
opaque

REF 414.10  
REF 414.11



## Cotton Wool Center

complete with insert with circular opening for cotton wool and narrow bars for small cotton pellets and pads.

smoked  
opaque

REF 414.40  
REF 414.41



smoked

## Cotton Wool Center

complete with insert with circular opening for cotton wool and broad bars for larger cotton pellets and pads.

smoked  
opaque

REF 414.50  
REF 414.51



## Cotton Roll Container

for the hygienic storage of cotton rolls in different sizes. The removable dividers permit numerous applications.

smoked  
opaque

REF 415  
REF 415.11



**!** Cotton wool, rolls and pellets – all perfectly packaged. Our high-quality systems are very versatile. They are made of solid acrylic glass and will **not lose their shape, even when heated up to 90 °C**. They are also extremely compact, measuring just 181 x 109 x 80 mm. An overlapping lid allows you to operate the dispenser with one hand. Our systems are also ideal for use in proximity to patients in accordance with the relevant RKI guideline (Germany).



⚠ Simply hygienic: our cotton wool dispensers are easy to disinfect with a surface disinfectant. We recommend our disinfectant sprays, Bechtozid plus/ Bechtozid sensitive (material compatibility proven by Röhm-Test).



## Container with lid

with bars for the universal use of cotton pellets and pads.

Size: 97 x 68 x 68 mm

smoked	<b>REF 417.10</b>
opaque	<b>REF 417.11</b>

## Container with lid

with circular opening for cotton wool and narrow bars for small cotton pellets and pads.

smoked	<b>REF 417.40</b>
opaque	<b>REF 417.41</b>



## Container with lid

with circular opening for cotton wool and broad bars for larger cotton pellets and pads.

smoked	<b>REF 417.50</b>
opaque	<b>REF 417.51</b>



⚠ Just right – for your individual requirements: Our cotton centers and cotton roll containers are available in opaque or smoked acrylic glass.

## Cotton Roll Container

with divider for two different sizes of cotton rolls.

smoked	<b>REF 418</b>
opaque	<b>REF 418.11</b>



## Derfla-Card-Plus

It's easy to handle and offers even more possibilities – and it helps you to document all important prophylactic procedures in detail: the Derfla-Card-Plus is the perfect assistant in your everyday work.

You can use the cards to note the Periodontal Screening Index (PSI) and prophylactic procedures for children.

The format and alphabetic tabs of the Derfla-Card-Plus are compatible with all standard card index systems. The way in which periodontal status and treatment schedule are recorded on the cards corresponds to the forms used by public health insurance funds (Germany).

The cards also have a heading for existing and planned dental prostheses and an easy-to-follow bookkeeping section on the back. Documents and X-ray films can be stored safely in the transparent sleeve provided.



- Optimum support: Offer young patients from 30 months to 18 years all the necessary prophylactic treatment covered by public health insurance.
- Prophylactic procedures covered by public health insurance are not subject to budgeting (only in Germany applies).

**Parodontalstatus und Behandlungsplan**

**Parodontaler Screening Index PSI**

**PSI-Erläuterung:**

Code	1	2	3	4	5	6
Grad	1-3,4	1-3,4	1-3,4	1-3,4	1-3,4	1-3,4
Art	gerade	gerade	gerade	gerade	gerade	gerade
Diagnose	gingivale	gingivale	gingivale	gingivale	gingivale	gingivale
	Parodontitis initial	Parodontitis initial	Parodontitis initial	Parodontitis initial	Parodontitis initial	Parodontitis initial

**Prophylaxe**

**Prophylaxekarte**

**Früherkennungsuntersuchung:** Vom 30. bis 72. Lebensmonat drei Mal abrechenbar, nach je 12 Monaten.

**Fluoridierung:** Vom 30. bis 72. Lebensmonat jedes Kalenderhalbjahres abrechenbar bei hohem Kariesrisiko (dmf-t)

**Individualprophylaxe:** Vom 6. bis 18. Lebensjahr abrechenbar.

**Mundhygienestatus:** (Vordruck) Abrechenbar jedes Kalenderhalbjahr, mit mindestens vier Monaten Zwischenraum. Bei 12- bis 19-jährigen ins Bonusheft.

**Mundgesundheitsaufklärung:** z.B. Putztechnik, Zahnselbst, Ernährung, Fluoridierung etc. Abrechenbar jedes Kalenderhalbjahr.

**Fluoridierung:** Abrechenbar jedes Kalenderhalbjahr. Bei hohem Kariesrisiko (dmf-t) zwei Mal pro Kalenderhalbjahr.

**Versiegelung der Zähne 6 und 7:** So oft wie klinisch indiziert.

**Hohes Kariesrisiko wenn:**

Alter bis	3 Jahre	dmf-t > 0
4 Jahre	dmf-t > 2	
5 Jahre	dmf-t > 4	
6 Jahre	dmf-t > 5	

Der dmf-t-Index ist ein internationaler Index zur Feststellung epidemiologischer Daten. Man nutzt ihn, um darzustellen, wie das Gebiss durch Zahnkrankheiten bislang beeinträchtigt wurde.

**d** = decayed = kariös  
**m** = missing = fehlend  
**f** = filled = gefüllt  
**t** = tooth = Zahn

Für das bleibende Gebiss werden Großbuchstaben (DMFT), für das Milchgebiss kleine Buchstaben (dmft) verwendet.

**Material-Übersicht**

Zahn	Produkt/Hersteller	Zahnfarbe/Zahnform
		Zahnfarbe OK
		Zahnfarbe UK
		Zahnform

Derfla-Card-Plus, complete with cover card, patient's record card, filing device, transparent envelope, file headings.  
Size: 209 x 155 mm

	System Card	Cover Card
grey	REF 60.1	REF 61.1
yellow	REF 60.2	REF 61.2
red	REF 60.3	REF 61.3
green	REF 60.4	REF 61.4
blue	REF 60.5	REF 61.5

Pack with 100 pieces



◀◀ **Patient's Record Card**

white, for the public health system  
white, for private insurances  
Pack with 100 pieces

**REF 72**  
**REF 72.1**



**Transparent Envelope**

Pack with 100 pieces

▲ **REF 73**

**File Headings**

color selection upon request

Pack with 100 pieces **REF 75**

**Filing Device**

self-adhesive

Pack with 100 pieces **REF 74**



**Record Cards**

They always find high acceptance. With tooth diagram and practical divisions printed on both sides for keeping clear accounts. Made of finest cardboard.

Size: 210 x 148 mm

- white
- yellow
- red
- green
- blue

- REF 89.1**
- REF 89.2**
- REF 89.3**
- REF 89.4**
- REF 89.5**

Pack with 100 pieces

**Transparent envelope**

(without illustration)

for record card, REF 89, sealed on three sides

Pack with 100 pieces

**REF 69**



## Mounts for X-Ray Films ▶▶

green cardboard, with cellophane windows for films size 2 x 3 cm or 3 x 4 cm. Protect and hold films safely.

Size: 14,8 x 10,5 cm

for **1 film**,

Pack with 100 pieces

**REF 81**

for **2 films**,

Pack with 100 pieces

**REF 82**

for **4 films**,

Pack with 50 pieces

**REF 84**

for **6 films**,





Pack with 50 pieces

**REF 86**

for **11 films**, size: 14,8 x 21 cm

Pack with 25 pieces

**REF 83.11**




## Mounts for X-Ray Films

envelope type, with white identification cards and cellophane window for protecting films, size 2 x 3 cm, 3 x 4 cm and 4 x 5 cm.

Size of card: 6,8 x 8,3 cm  
Pack with 100 pieces

REF 85



## Endless File for X-Ray Films

made of special plastic foil, transparent with white marking strip.

Max. size of films: 2 x 3 cm or 3 x 4 cm.

Pack with 160 pockets

REF 87

Pack with 480 pockets

REF 87.1

- ▼ Easy to separate the requested number of pockets
- ▼ Durable marking with all types of ball-point pens

## Panoramic Bags

of hard PVC sheet, transparent, neutral to films, for one film 13 x 18 cm, 18 x 24 cm, 12 x 30 cm or 15 x 30 cm, with 4 hole perforation for individual filing. May also receive a guide allowing them to be stored in a hanging cabinet.

Pack with 50 pieces

REF 88

## Self-adhesive Labels

with field for address, for panoramic bags.


Pack with 50 pieces

REF 88.1







**BEST**  
seller




white



blue



green



yellow

## Topliner Tray Paper

Our tray paper is made of highly absorbent crepe paper and chlorine free bleached. Available in different colors and sizes.

Size: 280 x 180 mm, dispenser box

white  
blue  
green  
yellow


REF 30.1  
REF 30.2  
REF 30.3  
REF 30.4

Size: 280 x 360 mm, blistered pack

white  
blue  
green  
yellow

REF 31.1  
REF 31.2  
REF 31.3  
REF 31.4

Contents: 250 pieces/pack

 The Topliner tray paper is available in four colors. Bring color into your dental office and choose your favorite.



The helpful assistant with plenty of bite: our articulating paper is colored on both sides and has white contact ends. The paper changes to a darker color to indicate the depth of bite. The paper has a thickness of 80 µ.



### Articulating Paper

straight, blue/blue  
straight, blue/red

Box with 12 books with 12 sheets

REF 407

REF 407.1

**BEST**  
seller



### Articulating Paper

Horseshoe type, blue/blue  
Horseshoe type, blue/red

Box with 12 books with 12 sheets

REF 405

REF 405.1



### Universal Matrix Strips

of acetate foil, 0,06 mm  
Size: 8 x 140 mm

Box of 100 pieces

REF 437

### Crown Setting Instrument


made of acrylic glass, not autoclavable, ergonomic shape, handy and solid type. For cleaning use a surface disinfectant.



REF 412



**BEST**  
seller



### Mixing Bowls

made of PVC, blue, easy to clean, suitable for plaster and impression compounds.

small, 100 mm in dia., 80 mm high,  
capacity: 300 ml

REF 454.1

large, 130 mm in dia., 100 mm high,  
capacity: 600 ml

REF 455.1



Do you find you acquiring ever more polishers, brushes and geometries – often with identical characteristics and the same composition? ProfiLab takes a new approach!  
 ProfiLab tidies things up: we offer carefully coordinated, tested and user-friendly product ranges to cover all standard materials and polishing levels. The benefit for you is that ProfiLab products achieve fast polishing results, durability and excellent performance in return of your investment.

## Polishing Paste

Universal paste for pre-polishing of precious alloys, acrylics and composites. This polishing paste can also be used as an alternative to pumice powder. Its soft consistency makes it easy to measure out and it causes less contamination of the working area than pumice products.

Block of 200 g

REF 192



## Polishing Emulsion

A high-shine emulsion for final polishing of precious alloys, acrylics and composites. Excellent results in a short time period. After polishing the material has reached a smooth and sealed surface. Best results in combination with Becht polishing wheels, REF 292.1. Fresh scent, skin friendly, non-toxic and economical.

Bottle of 100 ml

REF 191

## Pre-polishing of Ceramics/Composites/Zirconium



REF	230.1	230.2	230.4
Application	Multi-purpose polisher for ceramics and composites	Multi-purpose polisher for ceramics and composites	Multi-purpose polisher for ceramics and composites
RPM	4.000 - 10.000 min-1	4.000 - 10.000 min-1	4.000 - 10.000 min-1
Product description	Silicone polisher, wheel, unmounted	Silicone polisher, lens, unmounted	Silicone polisher, pointed, unmounted
Size	22 x 3 mm	22 x 4 mm	6 x 24 mm
Packing unit	Box of 12 pieces	Box of 12 pieces	Box of 12 pieces

## ProfiLab Brushes and Polishers



REF	249.1	249.2	249.3
Application	Polisher for ceramics and zirconium	Polisher for ceramics and zirconium	Polisher for ceramics and zirconium
RPM	4.000 - 10.000 min-1	4.000 - 10.000 min-1	4.000 - 10.000 min-1
Product description	Ceramically bonded polisher, cylinder, with diamond grains, mounted for HP	Ceramically bonded polisher, cylinder, with diamond grains, mounted for HP	Ceramically bonded polisher, cone, with diamond grains, mounted for HP
Size	6 mm Ø	4 mm Ø	4 mm Ø
Packing unit	piece	piece	piece



REF	1140.1	1141.1	1142.1
Application	Pre-polishing of ceramics	Pre-polishing of ceramics	Pre-polishing of ceramics
RPM	7.000 - 12.000 min-1	7.000 - 12.000 min-1	7.000 - 12.000 min-1
Product description	Diamond polisher, flame, coarse, mounted for HP	Diamond polisher, wheel, coarse, mounted for HP	Diamond polisher, disk, coarse, mounted for HP
Size	4 x 13 mm	11 x 2 mm	19 x 0,5 mm
Packing unit	piece	piece	piece

## Polishing of Ceramics



REF	1146.1	1147.1	1148.1
Application	Polishing of ceramics	Polishing of ceramics	Polishing of ceramics
RPM	7.000 - 12.000 min-1	7.000 - 12.000 min-1	7.000 - 12.000 min-1
Product description	Diamond polisher, flame, medium, mounted for HP	Diamond polisher, wheel, medium, mounted for HP	Diamond polisher, disk, medium, mounted for HP
Size	4 x 13 mm	11 x 2 mm	19 x 0,5 mm
Packing unit	piece	piece	piece

## High-shine polishing of Ceramics and Composites



REF	246.1	247.1
Application	High-shine polishing of ceramics	High-shine polishing of ceramics
RPM	7.000 - 12.000 min-1	7.000 - 12.000 min-1
Product description	Diamond polisher, flame, fine, mounted for HP	Diamond polisher, wheel, fine, mounted for HP
Size	4 x 13 mm	11 x 2 mm
Packing unit	piece	piece



REF	104	109.1	110.1
Application	High-shine polishing of ceramics	High-shine polishing of ceramics	High-shine polishing of ceramics
RPM	1.500 min-1	1.500 min-1	1.500 min-1
Product description	Polishing felt, disk, unmounted	Polishing felt, flame, mounted for HP	Polishing felt, pointed, mounted for HP
Size	21 x 3 mm	9 x 13 mm	7 x 15 mm
Polishing	High-shine Emulsion	High-shine Emulsion	High-shine Emulsion
Packing unit	Box of 12 pieces	Box of 6 pieces	Box of 6 pieces



REF	293.1	290.1	292.1	296.1
Application	High-shine polishing of composites	High-shine polishing of ceramics and composites	High-shine polishing of ceramics and composites	High-shine polishing of ceramics and composites
RPM	15.000 min-1	5.000 - 8.000 min-1	5.000 - 8.000 min-1	5.000 - 8.000 min-1
Product description	Goat hair, geared, mounted for HP	Linen, wheel, mounted for HP	Cotton, wheel, mounted for HP	Microfiber leather, wheel, mounted for HP
Size	22 mm Ø	22 mm Ø	22 mm Ø	22 mm Ø
Polishing	High-shine Emulsion/ProfiLab Polishing Paste	High-shine Emulsion	High-shine Emulsion	High-shine Emulsion
Packing unit	Box of 12 pieces	Box of 12 pieces	Box of 12 pieces	Box of 12 pieces



**BEST**  
seller

REF	296	279	1100.1 Taurus	1112.1 Taurus
Application	High-shine polishing of composites	High-shine polishing of composites	High-shine polishing of ceramics and composites	High-shine polishing of ceramics and composites
RPM	1.400 min-1	1.400 min-1	max. 16.000 min-1	max. 16.000 min-1
Product description	Microfiber polishing buffer, leather, 7 plies	Polishing buffer, soft tissue, 15 plies	Goat hair, wheel, mounted for HP	Horse hair, wheel, mounted for HP
Size	100 mm Ø	100 mm Ø	19 mm Ø	19 mm Ø
Polishing	High-shine Emulsion	High-shine Emulsion	High-shine Emulsion ProfiLab Polishing Paste	High-shine Emulsion ProfiLab Polishing Paste
Packing unit	piece	piece	Box of 12 pieces	Box of 12 pieces



REF 138.3

## Pre-polishing of Acrylics

**BEST**  
seller



REF	137/138	138.3	139
Application	Pre-polishing of acrylics	Pre-polishing of acrylics	Pre-polishing of acrylics
RPM	2.800 min-1	2.800 min-1	2.800 min-1
Product description	Brush wheel, horse hair/Brush wheel, Chungking bristles	Brush wheel, Chungking bristles, insert of muslin	Brush wheel, goat hair
Size	48 mm Ø	48 mm Ø	48 mm Ø
Polishing	Pumice powder/ProfiLab Polishing Paste	Pumice powder/ProfiLab Polishing Paste	Pumice powder/ProfiLab Polishing Paste
Packing unit	Box of 12 pieces	Box of 12 pieces	Box of 12 pieces

# ProfiLab Brushes and Polishers



REF	158.1	159.1	160.1	170.1 straight 171.1 convergent
Application	Pre-polishing of acrylics	Pre-polishing of acrylics	Pre-polishing of acrylics	Pre-polishing of acrylics
RPM	1.400 - 2.800 min-1	1.400 - 2.800 min-1	1.400 - 2.800 min-1	1.400 - 2.800 min-1
Product description	Brush wheel, Chungking bristles, 1 row	Brush wheel, Chungking bristles, 2 rows, convergent	Brush wheel, Chungking bristles, 3 rows, convergent	Brush wheel, Chungking bristles, 4 rows
Size	25/44 mm Ø	25/50 mm Ø	28/55 mm Ø	40/80 mm Ø
Polishing	Pumice powder/ProfiLab Polishing Paste	Pumice powder/ProfiLab Polishing Paste	Pumice powder/ProfiLab Polishing Paste	Pumice powder/ProfiLab Polishing Paste
Packing unit	piece	piece	piece	piece



REF 170.1



REF 171.1



REF 174.3

**BEST  
seller**



REF	173.3	174.3	176/176.1
Application	Pre-polishing of acrylics	Pre-polishing of acrylics	Pre-polishing of acrylics
RPM	1.400 min-1	1.400 min-1	1.400 - 2.800 min-1
Product description	Brush wheel, Chungking bristles and muslin insert	Brush wheel, Chungking bristles and leather insert	Brush wheel, white nylon bristles and fleece insert
Size	40/80 mm Ø	40/80 mm Ø	25/55 mm Ø / 40/80 mm Ø
Polishing	Pumice powder/ProfiLab Polishing Paste	Pumice powder/ProfiLab Polishing Paste	Pumice powder/ProfiLab Polishing Paste
Packing unit	piece	piece	piece



REF	1130.1	1132.1	1133.1
Application	Pre-polishing of acrylics	Pre-polishing of acrylics	Pre-polishing of acrylics
RPM	1.400 - 10.000 min-1	1.400 - 10.000 min-1	1.400 - 10.000 min-1
Product description	Scotch Brite™, Aluminiumoxyde, coarse, mounted for HP	Scotch Brite™, Aluminiumoxyde, fine, mounted for HP	Scotch Brite™, Aluminiumoxyde, super fine, mounted for HP
Size	25 mm Ø	25 mm Ø	25 mm Ø
Packing unit	Box of 6 pieces	Box of 6 pieces	Box of 6 pieces



**BEST  
seller**

REF	294.1	1100.1 Taurus	1112.1 Taurus
Application	Pre-polishing of acrylics	Pre-polishing of acrylics	Pre-polishing of acrylics
RPM	15.000 min-1	max. 16.000 min-1	max. 16.000 min-1
Product description	Bristles, geared, mounted for HP	Goat hair, wheel, mounted for HP	Horse hair, wheel, mounted for HP
Size	22 mm Ø	19 mm Ø	19 mm Ø
Polishing	Pumice powder/ProfiLab Polishing Paste	Pumice powder/ProfiLab Polishing Paste	Pumice powder/ProfiLab Polishing Paste
Packing unit	Box of 12 pieces	Box of 12 pieces	Box of 12 pieces




REF	236.2	235.2	1237.2
Application	Pre-Polishing of acrylics	Pre-Polishing of acrylics	Pre-Polishing of acrylics
RPM	2.000 - 7.000 min-1	2.000 - 7.000 min-1	2.000 - 7.000 min-1
Product description	Silicone polisher, cone, medium, mounted for HP	Silicone polisher, cone, coarse, mounted for HP	Silicone polisher, cone, fine, mounted for HP
Size	10 x 24 mm	10 x 24 mm	10 x 24 mm
Packing unit	Box of 6 pieces	Box of 6 pieces	Box of 6 pieces

## High-shine polishing of Acrylics



REF 299



REF	279	296	299
Application	High-shine polishing of acrylics	High-shine polishing of acrylics	High-shine polishing of acrylics
RPM	1.400 min-1	1.400 min-1	1.400 min-1
Product description	Polishing buffer, soft tissue, 15 plies	Microfiber polishing buffer, leather, 7 plies	Polishing buffer, muslin, 30 plies
Size	100 mm Ø	100 mm Ø	100 mm Ø
Polishing	High-shine Emulsion	High-shine Emulsion	High-shine Emulsion
Packing unit	piece	piece	piece



REF	290.1	292.1
Application	High-shine polishing of acrylics	High-shine polishing of acrylics
RPM	5.000 - 8.000 min-1	5.000 - 8.000 min-1
Product description	Linen, wheel, mounted for HP	Cotton, wheel, mounted for HP
Size	22 mm Ø	22 mm Ø
Polishing	High-shine Emulsion	High-shine Emulsion
Packing unit	Box of 12 pieces	Box of 12 pieces



REF	453.1	453.2	453.3
Application	Mixing of acrylics	Mixing of acrylics	Mixing of acrylics
Material	of silicone	of silicone	of silicone
Product description	Small bowl, 32 mm ø, 28 mm high, suction cup, contents: 5 ml, set of 3 pieces	Medium bowl, 50 mm ø, 45 mm high, suction cup, contents: 25 ml	Large bowl, 70 mm ø, 60 mm high, suction cup, contents: 75 ml

# ProfiLab Brushes and Polishers

## Pre-polishing of Alloys and Castings



REF	230.1	230.2	230.4
Application	Pre-polishing of precious alloys and castings	Pre-polishing of alloys and castings	Pre-polishing of alloys and castings
RPM	4.000 - 10.000 min-1	4.000 - 10.000 min-1	4.000 - 10.000 min-1
Product description	Silicone polisher, wheel, coarse, unmounted	Silicone polisher, lens, coarse, unmounted	Silicone polisher, pointed, coarse, unmounted
Size	22 x 3 mm	22 x 4 mm	6 x 24 mm
Packing unit	Box of 12 pieces	Box of 12 pieces	Box of 12 pieces



REF	237.1	237.2	237.3
Application	Pre-polishing of alloys and castings	Pre-polishing of alloys and castings	Pre-polishing of alloys and castings
RPM	7.000 - 12.000 min-1	7.000 - 12.000 min-1	7.000 - 12.000 min-1
Product description	Silicone polisher, wheel, coarse, unmounted	Silicone polisher, lens, coarse, unmounted	Silicone polisher, pointed, coarse, mounted for HP
Size	22 x 3 mm	22 x 4 mm	5 x 16 mm
Packing unit	Box of 12 pieces	Box of 12 pieces	Box of 12 pieces

## Polishing of Alloys and Castings



**BEST  
seller**

REF	137 Horse Hair 138 Chungking-Bristles	178.1	179.1
Application	Polishing of alloys and castings	Polishing of alloys and castings	Polishing of alloys and castings
RPM	2.800 min-1	2.800 min-1	2.800 min-1
Product description	Bristle brush	Brush wheel, short Chungking bristles, 2 rows, convergent	Brush wheel, short Chungking bristles, 2 rows, convergent
Size	48 mm Ø	40/60 mm Ø	50/70 mm Ø
Polishing	ProfiLab Polishing Paste	ProfiLab Polishing Paste	ProfiLab Polishing Paste
Packing unit	Box of 12 pieces	piece	piece



**BEST  
seller**

REF	294.1	1100.1 Taurus	1112.1 Taurus
Application	Polishing of alloys and castings	Polishing of alloys and castings	Polishing of alloys and castings
RPM	15.000 min-1	max. 16.000 min-1	max. 16.000 min-1
Product description	Bristle brush, geared, mounted for HP	Goat hair, wheel, mounted for HP	Horse hair, wheel, mounted for HP
Size	22 mm Ø	19 mm Ø	19 mm Ø
Polishing	ProfiLab Polishing Paste	ProfiLab Polishing Paste	ProfiLab Polishing Paste
Packing unit	Box of 12 pieces	Box of 12 pieces	Box of 12 pieces

## High-shine polishing of Alloys and Castings



REF	290.1	292.1	296.1	1125
Application	High-shine polishing of alloys and castings	High-shine polishing of alloys and castings	High-shine polishing of alloys and castings	High-shine polishing of alloys and castings
RPM	5.000 – 8.000 min-1	5.000 – 10.000 min-1	5.000 – 8.000 min-1	5.000 – 10.000 min-1
Product description	Linen, wheel, mounted for HP	Cotton, wheel, mounted for HP	Microfiber leather, wheel, mounted for HP	Polishing disc, cardboard, coated, <b>unmounted</b>
Size	22 mm Ø	22 mm Ø	22 mm Ø	25 mm Ø
Polishing	High-shine Emulsion	High-shine Emulsion	High-shine Emulsion	ProfiLab Polishing Paste
Packing unit	Box of 12 pieces	Box of 12 pieces	Box of 12 pieces	Pack of 100 pieces



REF 299



REF	279	296	299
Application	High-shine polishing of alloys and castings	High-shine polishing of alloys and castings	High-shine polishing of alloys and castings
RPM	1.400 min-1	1.400 min-1	1.400 min-1
Product description	Polishing buffer, soft tissue, 15 plies	Microfiber polishing buffer, leather, 7 plies	Polishing buffer, muslin, 30 plies
Size	100 mm Ø	100 mm Ø	100 mm Ø
Polishing	High-shine Emulsion	High-shine Emulsion	High-shine Emulsion
Packing unit	piece	piece	piece



REF	1704	1705
Product description	Casseroles, of porcelain	Casseroles, of porcelain
Application	Removal of oxide layers	Removal of oxide layers
Size	75 mm Ø	85 mm Ø

# Index

	Page		Page		Page			
A	AirFresh med	21	Diamond Cleansing Brush	22	ProfiDent	12, 13		
	Alcohol lamps	35, 36	Disinfection and Cleaning Liquids, Sundries	14 – 22	ProfiGuard/ProfiPolish	8		
	AquaHy	21			ProfiLab	47 – 53		
	Articulating Paper	45			Prophy Cups, latex-free	11		
	Asphalin B Tablets	21	E	Endodontics	4 – 7	Prophylaxis	8 – 11	
	Asphalin Trays	33	F	Felts	49	Prophylaxis brushes	10, 11	
	Asphaline temp	4		Film Cards	42	ProphyRepair	9	
				Fräsator	30			
B	Becht Carborundol	16	G	Glass Blocks and Slabs	37	R	Record Cards	40, 41
	Bechtid plus	17		Guttapercha Points	7		Rinsing Jar, crystal glass	36
	Bechtifix alcohol-free	20	I				Rinsing Jar, ÖKOPAL	36
	Bechtifix plus	20		Instrument Rest „Antifekt“	27	S	Scotch Brite® Brushes	50
	Bechtol	16		Instrument Trays	27		Silicone Membranes	23
	Bechtol CLASSIC	15	L				Standard Trays	27
	Bechtol plus	15		Laboratory brushes for acrylics	49 – 51		Stainless Steel Lids	29
	Bechtol Premium	14		Laboratory brushes for ceramics, composites, zirconium, alloys	47 – 49		Steel Separating Strips	12
	Bechtoid plus	18		Lids, stainless	29	T	TAURUS Brushes	49, 50, 52
	Bechtoid Premium	18	M				Tooth Cleaning Brushes	11
	Bechtoid sensitive	19		Matrix Strips	46		Topliner Tray Paper	44
	Bechtoid soft	17		Medicament Bottles	24, 25		Tweezers	30, 33
	Bottle Stands	26		Medi-Spencers	23	W	Waste Receiver	29, 31
	Bur Cleansing Brushes	22		Mixing Bowls	46, 51	X	X-Ray Record Systems	43
			Mounds for X-Ray Films	42, 43				
C	Casserole	53	O	Office Administration	40 – 43			
	Closure UNIKAP	28		Opaline and Glass Products	23 – 37			
	Combi Stand	23	P					
	Cotton Center, acrylic	38		Panoramic Bags	43			
	Cotton Container, acrylic	39		Paper Points	6			
	Cotton Container, glass	31, 32		Paper Tissues	26			
	Cotton Jars	32		Paromatic-brush	10			
	Crown Setting Instrument	46		Pellet Container	31, 32			
	Cup Brushes, prophylaxis	11		Perma Evolution	5			
				Petri Dishes	33, 34			
D	Dappen Dishes	24		Polishing Emulsion	47			
	Dental Floss	10		Polishing Paste	47			
	Dental Polishers, chair side	12 – 13		Polishing Strips, Amalgam	12			
	Denture Bowl	28						
	Derfla-Card-Plus	40						
	Disinfecting Jar	29						
	Disinfection Tub	22						



## Contact Persons

**Ulrich A. Becht**

Managing Director

Tel. +49 781 60586-0

becht@becht-online.de

**Harald Kolloch**

Vice-President, Sales and Marketing

Tel. +49 781 60586-0

kolloch@becht-online.de

**Susanne Buehler**

Assistant Sales and Marketing Export

Tel. +49 781 60586-21

buehler@becht-online.de

**Katrin Stokert**

Assistant Sales Export

Tel. +49 781 60586-24

stokert@becht-online.de

**Elvira Meier**

Marketing

Tel. +49 781 60586-12

meier@becht-online.de

**Christian Gottschling**

Area Export Manager

Tel. +49 331 23700347

Mobile: +49 151 55569900

gottschling@becht-online.de

**Linda Yang 杨琳丹**

Market Development Manager CHINA

中国区市场拓展经理

Mobile: +86 18600500478

linda@becht-online.de

**Angelika Noel**

Quality Management Systems

Tel. +49 781 60586-17

noel@becht-online.de

**Bettina Glatz**

Accounting

Tel. +49 781 60586-30

glatz@becht-online.de



ALFRED BECHT GMBH

BESSER BECHT. MADE IN OFFENBURG.



Alfred Becht GmbH  
P.O.B. 11 45  
D-77601 Offenburg

Tel. +49 781 60586-0  
Fax +49 781 60586-40  
[www.becht-online.de](http://www.becht-online.de)



ALFRED BECHT GMBH

BESSER BECHT. MADE IN OFFENBURG.