

All you need for easy lab work



SpofaDental
A Kerr Company

Partnering your practice.

crown & bridge

Superpont C+B
Superpont
Duracryl
Conalor
Spofacryl



Superpont C+B

Crown and bridge acrylic resin

- direct modelling of acrylic veneered crowns and bridges
- acrylic jacket crowns
- temporary crowns and bridges

high modelling capability



Superpont

Heat-curing crown and bridge acrylic resin

- acrylic jacket crowns and bridges
- temporary crowns and bridges

long-term proven colour stability

Duracryl

Self-curing crown resin

- manufacturing of temporary crowns, bridges
- repairing acrylic crowns and veneered crowns

short polymerization time



Conalor

Full-scale system of universal opaque colours for dyeing the metal frameworks of acrylic veneered fixed restorations

- dyeing under denture base resin on saddles of framework of removable partial dentures
- toning and shading of crown resins
- imitation of enamel cosmetic defects



very good metal frameworks coverings

Spofacryl

Self-curing adhesive resin

- affixing loosed veneers in fixed restorations
- repairs of involved occlusal surfaces, temporary splinting of traumatically damaged teeth
- affixing fixed orthodontic appliances



excellent plasticity, easy slumping



Casting

Silikan
Silisan N
Gloria Special
KP50



Silikan

Phosphate-bonded investment material

- high-melting alloys
- manufacturing of cast models
- manufacturing of cast moulds for the metal frameworks of partial removable dentures and bar constructions.

high mould strength



Silisan N

Expansion mixing solution for Silikan

- ensuring optimal properties of Investment material

optimization of phosphate-bonded investment material properties

Gloria Special

Gypsum-bonded investment material

- for low-melting dental alloys

very good reproduction and optimum dimensional change



KP50

Ceramic casting ring liner

- enabling total expansion of investment materials

exact casting preventing formation of cracks, asbestos-free

Prosthesis making

Tessex
Duracrol
Ceradent
Stomaflex Varnish
Dentiplast
Teeth
Superacryl Plus
Duracryl Plus



Tessex

Shellac base plates for individual impression trays

- manufacturing impressions of edentulous jaws
- manufacturing of wax bite registration templates

also with aluminium admixture



Duracrol

Self-curing acrylic resin

- fabrication of individual impression trays
- manufacturing of individual phantom models

easy handling and finishing

Ceradent

Base plate modelling wax

- broad range of use in dental laboratory.

easy formable



Stomaflex Varnish

Silicone-based separating material

- creates necessary separating and affixative layer of interproximate spaces of teeth prior to manufacturing of removable restoration by flasking.

ideal application, no slumping

Dentiplast

Auxiliary impression material from thermoplastic wax-resin mass

- edge contouring of individual impression trays in mucodynamic impressions
- edge contouring of dentures prior to impressions for indirect rebasing
- mucodynamic impressions to determine neutral zone of the complete dentures

great adapting



Teeth

Three-layer, cross-linked acrylic teeth

- unique design satisfies all demands of modern dentistry on artificial teeth
- established with sense of maximum aesthetic details, comfort and above all functionality
- characterized by colour stability, high resistance to abrasion and perfect junction with base resin
- delivered in VITA™ shades and aligned according to innovated combination of occlusal relationships.

long-lasting colour stability



Superacryl Plus

Heat-curing denture base resin

- total dentures
- partial dentures

durable adhesion to acrylic teeth.



Duracryl Plus

Self-curing denture base resin

- repair of removable dentures
- manufacturing of immediate dentures

durable bonding to denture base resin



Orthodontic

Premacryl Plus



Premacryl Plus

Self-curing orthodontic acrylic resin

- manufacturing of orthodontic appliances
- immediate dentures
- repairs of removable dentures.

wide range of colour shades



Stones

Convertin® Hart
Isodent
Mramorit
Alabaster Plaster
Efektor



Convertin® Hart

First-class quality dental stone – type 4

- manufacturing of detail situation models in fixed prosthetics
- high-strength gypsum conforming demands on reproduction of the finest details with high flowability.

low expansion; high marginal strength



Isodent

Gypsum separating solution

- provides perfect isolation of gypsum models and moulds
– gypsum/gypsum, gypsum/resin

easy application



Mramorit

High quality, proven model dental stone
– type 3

- preparation of moulds for removable restorations,
- fixed restorations
- fixation of models into articulators

available also in white



Alabaster Plaster

Model plaster – type 2

- suitable for manufacturing of situation models
- moulds for removable and fixed restorations
- affixation of models into articulators.

wide range of use in dental laboratory



Efaktor

Fast setting model plaster – type 2

- affixation of models into articulators

satisfying all demands on auxiliary works
in dental laboratory

Order information

Making of crown & bridge

Superpont® C+B Box 4324150

8 x 20 g of cervical shades: A3.5, A4, B2, B3, C4, D2, 2xD4; 14 x 20 g of dentin shades: A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C3, C4, D2, D3, D4; 7 x 20 g of enamel shades: A2, 2xA3.5, B2, C4, D2, D4; 1 x 20 g of transparent shade: T; 2 x 50 g of liquid

Refills

Superpont C+B 100 g of powder

powder 100 g cervical A3,5	4324212
powder 100 g cervical B2	4324213
powder 100 g cervical B3	4324214
powder 100 g cervical A4	4324215
powder 100 g cervical C4	4324216
powder 100 g cervical D2	4324218
powder 100 g cervical D4	4324219
powder 100 g dentin A1	4324220
powder 100 g dentin A2	4324221
powder 100 g dentin A3,5	4324222
powder 100 g dentin B2	4324223
powder 100 g dentin B3	4324224
powder 100 g dentin A3	4324225
powder 100 g dentin C2	4324226
powder 100 g dentin C3	4324227
powder 100 g dentin D3	4324228
powder 100 g dentin D4	4324229
powder 100 g enamel A2	4324231
powder 100 g enamel A3,5	4324232
powder 100 g enamel B2	4324233
powder 100 g enamel C4	4324235
powder 100 g enamel D2	4324237
powder 100 g enamel D4	4324239
powder 100 g enamel transparent - T	4324240
powder 100 g dentin A4	4324261
powder 100 g dentin B4	4324262
powder 100 g dentin C4	4324263
powder 100 g dentin D2	4324265
250 g liquid VITA™	4324903

Superpont® Set 4321150

8 x 20 g of cervical shades: A3.5, A4, B2, B3, C4, D2, 2xD4; 14 x 20 g of dentin shades: A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C3, C4, D2, D3, D4; 7 x 20 g of enamel shades: A2, 2xA3.5, B2, C3, D2, D4; 1 x 20 g of transparent shade: T; 2 x 50 g of liquid

Refills

Superpont 100 g powder

powder 100 g cervical A3,5	4321212
powder 100 g cervical A4	4321213
powder 100 g cervical B2	4321215
powder 100 g cervical B3	4321216
powder 100 g cervical C4	4321217
powder 100 g cervical D2	4321218
powder 100 g cervical D4	4321219
powder 100 g dentin A1	4321221
powder 100 g dentin A2	4321222
powder 100 g dentin A3	4321223
powder 100 g dentin A3,5	4321224
powder 100 g dentin A4	4321225
powder 100 g dentin B2	4321226
powder 100 g dentin B3	4321227
powder 100 g dentin B4	4321228
powder 100 g dentin C2	4321232
powder 100 g dentin C3	4321233
powder 100 g dentin C4	4321234
powder 100 g dentin D2	4321236
powder 100 g dentin D3	4321237
powder 100 g dentin D4	4321238
powder 100 g enamel A2	4321241
powder 100 g enamel A3,5	4321242
powder 100 g enamel B2	4321244
powder 100 g enamel C4	4321246
powder 100 g enamel D2	4321248
powder 100 g enamel D4	4321249
powder 100 g transparent - T	4321251
250 g liquid VITA™	4321903

Duracryl® 4811100

7 x 10 g of powder in colour shades 1, 5, 9, 37, 45, V (charcoal grey), VI (brown)
50 ml of liquid

Conalor® 4813210

3 x 6 g of powder in VITA™ shades A, B, C (D) 4 x 3 g of powder in additional shades (white, brown, pink, grey)
16 g of liquid

Refills

15 g of powder, white	4813222
15 g of powder, brown	4813223
15 g of powder, pink	4813225
15 g of powder, grey	4813226
15 g of powder, A	4813231
15 g of powder, B	4813232
15 g of powder, C+D	4813233
25 g of liquid	4813250

Spofacryl®

100 g powder 0 + 100 g liquid	4318200
100 g powder VITA™ A2 + 100 g liquid	4318212
100 g powder VITA™ A3,5 + 100 g liquid 4	318213
100 g powder VITA™ B3 + 100 g liquid	4318223
100 g powder VITA™ C2 + 100 g liquid	4318232

Order information

Casting

Silikan

5 kg	4521120
10 kg	4521130

Silisan® N

1000 g	4531110
--------	---------

Gloria® Special

10 kg	4512110
-------	---------

KP 50

(Ceramic Casting Ring Liner)	4815210
20 m of tape 50mm wide	

Stones

Convertin Hart

1000 g sachet	4511210
5 kg bucket	4511215

Mramorit®

25 kg paper bag	4511135
-----------------	---------

Mramorit® Blue

25 kg paper bag	4511125
-----------------	---------

Alabaster Plaster

25kg paper bag	4511325
----------------	---------

Efektor®

25 kg paper bag	4251135
-----------------	---------

Isodent®

250 g of solution	4817220
-------------------	---------

Prosthesis making

Tessex

12 plates for upper jaw	4812112
12 plates for lower jaw	4812113

Tessex AL

12 plates for upper jaw	4812212
12 plates for lower jaw	4812213

Duracrol

powder 500 g	4331302
liquid 250 g	4331902

Ceradent®

1000 g of wax sheets, size 170 x 90 mm (aprox. 55 to 57 pcs)	4411115
---	---------

Stomaflex® Varnish

140 g of paste, 11,8 g of vulcanite	4817330
-------------------------------------	---------

Dentiplast®

20 g of paste	4232110
---------------	---------

Spofadent® Plus

Complete set 1 x 28 pcs	See order form
Posterior teeth 4 x 8 pcs	See order form
Anterior teeth 6 x 6 pcs	See order form
Complete set 2 x 28 pcs	See order form

Supercryl® Plus

500 g of powder O	4328411
500 g of powder U	4328412
500 g of powder V	4328413
500 g of powder Z	4328414
500 g of powder X	4328417
250 g of liquid	4328902

Duracryl® Plus

500 g of powder O	4316411
500 g of powder U	4316412
500 g of powder V	4316413
500 g of powder Z	4316416
250 g of liquid	4316902

Orthodontic

Premacryl Plus

powder 500 g O	4342400
powder 500 g V	4342405
250 g liquid O	4342921
250 g liquid R (red)	4342922
250 g liquid B (blue)	4342923
250 g liquid GR (green)	342924
250 g liquid G (yellow)	4342927

Our Dealers

SpofaDental
A Kerr Company

Partnering your practice.

SpofaDental a.s.

Markova 238, 506 01 Jicin, Czech Republic

tel.: +420 493 583 277, fax: +420 493 583 275

info@spofadental.com ■ www.spofadental.com